

# INFANT CARE TROLLEY



## Infant Trolley

- Easy access of infant from all sides with thick Acrylic collapsible side panels
- Acrylic baby tray with Infant's Head up / down facility & foam mattress
- Provision for positioning X-ray cassette
- S.S. Tray with slot to use undersurface Phototherapy
- I.V. Pole facility
- One lower shelf
- M.S. Fabrication with powder coated paint
- Mounted on castors for easy mobility
- Can be used along with Warmer stand, Phototherapy stand, Wall mounted Warmer & Phototherapy and as a separate bed for neonates
- Unit also available in Stainless Steel

Specification & design are subject to change without prior notice

9,500=00 GST EXTRA

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

6

© Copyright Reserved by : Western Surgical - Rajkot.