

Standard Model Intensive Care Warmer

SILENT FEATURES

EASY TO USE

As technology advanced, we have kept the Warmer simple & easy to use with one touch controls.

FOCUSED HEATING

The heater with reflector makes it easier to work with the baby by efficiently providing even heat distribution across the cradle and not into the surrounding environment.

AIR CONTROL MODE

Air control is an electronic feedback system which functions as a thermostat to maintain a constant temperature of baby cradle by regulating the heat output of Radiant Warmer. It Also acts as a WATCH DOG as Air over temperature cuts off heater event if the Warmer is operating in Servo mode.

MANUAL MODE

The Manual mode allows the Warmer to be kept 'baby ready' without nuisance alarms. When needed the Warmer provides safe, powerful & controlled warming to restore the temperature drop experienced by the baby immediately after delivery.

TECHNICAL SPECIFICATION

- Power Source : 230V AC, $\pm 10\%$, 50Hz.
- Light Source : LED Light to examine babies.
- Heater Capacity : 600 Watts.
- Power Consumption : 800 Watts at 100% Heater power.

SERVO CONTROL MODE

Servo control is an electronic feedback system which functions as a thermostat to maintain a constant temperature at the site of Skin Probe by regulating the heat output of Radiant Warmer maintaining a constant abdominal Skin temperature between 36.5 & 37.0° C is the simplest way to provide a "thermoneutral" environment, minimizing the no. of calories needed to maintain normal body temperature & reducing the risks of cold stress or overheating. Servo control is the only acceptable method of heat regulation for the infant cared for under a Radiant Warmer.

INFANT WARMER (SKIN / MANUAL)

Skin Temperature displayed on 1" bright Red LED display.
Set temperature displayed on separate bright Red LED display,
Two modes of Warming i.e. Skin, Air & Manual.
Heater output in 10 equal steps in Percentage.
Programmable Heater power & Time in Manual mode.
Individual LED for selection for mode (Skin, Air & Manual).

WSW-01



APGAR TIMER / ATTEND

BABY ALARM PROGRAMMABLE

MUTE TIME SET TEMPERATURE

RANG : SKIN 25° C-39° C

30,000=00 GST EXTRA

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

(2)

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.



Scanned with OKEN Scanner