

# WS-9001-T Touch Screen Multi-para Patient Monitor 12.1"



## FEATURE:

5 standard parameters: ECG, RESP, NIBP, SPO2, 1-TEMP.  
Real time S-T segment analysis, pace-maker detection and ARR analysis.  
Multi-display selectable including standard, large font, trend coaxis, OxyCRG dynamic  
On-line help and Patient info input management  
Multi-lead ECG waveforms display in phase. And capture dynamic waveforms.  
Large volume of tabular and graphic trends information storage easy to recall  
Efficient resistance to interference of defibrillator and HF knife.  
Up to 4 hours working capacity of built-in rechargeable battery.

## SPECIFICATIONS:

ECG: Lead mode: 5 -lead(R,L,F,N,C)  
Waveform: 3 and 7 channel selectable  
Heart rate range: Adult: 15~300bpm; Neonate:/pediatric:15~350bpm  
Resolution: 1bpm Filter: surgery mode:1~20Hz  
Alarm range: 15~1350bpm S-T segment detecton: Measurement range:2.0mV~+2.0mV  
Alarm Range: -2.0mV~+2.0mV  
Accuracy: -0.8mV ~+0.8Mv Error: +0.02Mv  
Arrhythmia analysis:YES

## SPO2:

Measurement range: 0-100%  
Resolution: 1%  
Accuracy: +2%(70-100%);0-69% unspecified  
Alarm range:0~100%  
Pluse Rate: range:20~300bpm  
Resolution: 1bpm  
Error:+1bpm or +2%/whichever is greater

## NIBP:

Method: Digital Automatic oscillometric  
Operation mode : Manual/Automatic/continuous  
Over~pressure protection: Double safety protection  
Resolution: 1 mmHg  
Alarm: Systolic.Diastolic,Mean

## RESPIRATION:

Method: Thoracic impedance  
Measurement range: Adult:7~120rpm;  
Nonate/Pediatric: 7~150rpm  
Apnea alarm: YES, 10~40's  
Resolution: 1rpm  
Accuracy: +2rpm

## TEMPERATURE:

Compatible probe: YES or CYF  
Measurement range: 5~50°C  
Accuracy: +0.1°C  
Refreshing time: about 1  
Average measuring time:<10s

Optional-Recorder, 8-IBP, 8-TEMP, 2-SpO2, Nellcor/ Masimo SpO2, Phaisein/ Resplronics Multi-gas.

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

(25)

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.



Scanned with OKEN Scanner