

Standard Model Power LED Phototherapy



- ▶ Time Totalizer (Optional).
- ▶ Power LED's emitting blue light in 450-470nm.
- ▶ Long life >45,000 hrs.
- ▶ Intensity : >25 μ w / cm² / nm.
- ▶ M.S. Fabrication with powder coated paint.
- ▶ Operation Voltage : 100-250 V AC, 50Hz.

12Pc. Hi power LED provides to cover the bed area and for effective Phototherapy treatment. No ultraviolet ray ensures minimize the UV exposures. No infrared rays ensures minimize the water loss. Adjustable Height - Conveniently allows you to control Radiance levels on the neonate and different types of cradle, Convenient Tilting of Source Unit, Height Adjustment with Single hand and Castor with elegant maneuvering helps to place the system at different positions, at user convenience. Lamp usages hours greater than 45,000 hours ensures no need of frequent change of lamp.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.