

INTENSIVE CARE INCUBATOR

Infant Incubator



MECHANICAL SPECIFICATION :

- ◆ Acrylic Canopy with front loading & 4 port Holes.
- ◆ Removable sub assemblies like Baby Tray, Blower System.
- ◆ Inlet for Oxygen, I.V. Tubing, Probes etc.
- ◆ Excellent access and sealing through gloves on port holes.
- ◆ Humidity tray & Humidity Indication by Hygrometer.
- ◆ Head Up / Down positioning facility & Foam Mattress.
- ◆ IV Stand Facility.
- ◆ Mounted on heavy duty castors for easy mobility.
- ◆ Fabrication in Mild Steel with Epoxy paint.

TECHNICAL SPECIFICATION :

- ◆ Highly reliable Micro Computer based temp. controller.
- ◆ Two modes of warming i.e. (Skin & Air).
- ◆ Individual LED for selection of mode (Skin & Air).
- ◆ Skin & Air Temperature displayed individually on 1" bright Red LED display.
- ◆ Set temperature displayed on separate bright Green LED display.
- ◆ Separate display indication for % of heater output in 10 equal steps
- ◆ Apgar Timer / Attend Baby Alarm / Programmable Mute Time.

SET TEMPERATURE RANGE :

Skin 25° C — 39° C Resolution : 0.1° C
Air 25° C — 39° C Accuracy : 0.1° C

SAFETY ALARMS (with Audio Visual message and Mute Facility)

Skin Mode :

High Temp. (> 1° C of Set temperature) HI
Low Temp. (< 1° C of Set temperature) LO
Over Temp (> 38° C) Ovr
SKIN PROBE FAIL : Prb

Air Mode :

High Temp. (> 1.5° C of Set temperature) HI
Low Temp. (- 3° C of Set temperature) LO
Over Temp (> 39° C) Ovr
AIR PROBE FAIL : Prb

Power Source : 230 V AC +/- 10%, 50 Hz.
Power Consumption : 600 Watts at 100% Heater power.
Heater Capacity : 400 Watts.
Power Fail : Audio Visual Alarm.

Model Available : with 3 drawers & Storage Cupboard.

Specification & design are subject to change without prior notice.

55,000=00 GST EXTRA

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.