


Contact Phototherapy



- ▶ Compact, Portable & Easy to Carry. ▶ Detachable baby Tray for Cleaning.
- ▶ No need to Cover eyes of the baby. ▶ Low Power Consumption (24W).
- ▶ LED's emitting blue light in 450-470 nm. Spectrum matching peak absorption wavelength (458nm.) at which bilirubin is broken down.
- ▶ Do not emit light in infrared radiation range-reducing the risk of fluid loss.
- ▶ Long Life > 20,000 hrs. ▶ Irradiance : > 50 μ w / cm² / nm.
- ▶ Can be used separately in mother bed / in trolley.
- ▶ M.S. Fabrication with powder coated paint.
- ▶ Operation Voltage : 100-250 V AC, 50Hz.
- ▶ Digital Time Totalizer for LED Usage & Therapy Time Resettable (Optional).

22,000=00 GST EXTRA

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

5

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.