



WS - M4



**FIXED HEIGHT ICU
CLASSIC MOTORIZED**

Size: 2100 x 965 x 460mm

- Motorized with 2 Pcs actuators 1 Control box,
1 handset with battery backup
- Backrest & Footrest Motorized.
- Trendelenberg & Reverse Trendelenberg on gas Spring
- Wheels in high quality engineering plastic
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high
strength engineering plastic
- 4 pcs uring hooks.
- 4 pole Saline attachment.

1,35,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

2

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.