



**WS - M1**



**Size 2150 x 965 x 510 - 730mm ht.**

### **ICU EXCEL MOTORIZED**

- Completely motorized 4 Pcs actuators 1 Control box, 1 Pc handset with battery backup
- Bed board in High strength engineering plastic
- Wheels in High quality engineering plastic With central Locking
- Head and foot end bows in high strength engineering plastic
  
- 4 underside occidental style bed side rails
- 4pcs urine hooks.
- 4 pole Saline attachment.
- Optionally available with urgent release (CPR) function
- Also Available with ACP if required

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**

1

[www.westernsurgical.in](http://www.westernsurgical.in)

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. : +91-281-2232871 / 2243159

© Copyright Reserved by : Western Surgical - Rajkot.