

Delivery Table Hydraulic



90,000/-

1. Overall Size : 1900mm L x 900mm W x 600-800mm H
2. Tabular Frame Work mounted on 125mm dia castor two with brakes.
3. Three Section top made of M.S Pipe with fixed cushioned mattress.
4. Trendlenburg / Reverse Trendlenburg Positions adjustment through Gear Mechanism.
5. Height Adjustment through Imported Hydraulic Pump.
6. Back section Adjustment through ratchet system.
7. Leg end section can slide under the main section.
8. Provided with a pair of Upholstery Lithotomy Positions
9. S.S Telescopic I.V Rod.
10. Detachable S.S Tray.



Electric

1. Overall Size : 1880mm L x 900 mm W x 600-800mm H
2. Tabular Frame Work mounted on 125mm dia castor two with brakes.
3. Three Section top made of M.S Pipe with fixed cushioned mattress.
4. Trendlenburg / Reverse Trendlenburg Positions , Height Adjustment, Backsection Adjustment through electric Actuator System.
5. Leg end section can slide under the main section.
6. Provided with a pair of Upholstery Lithotomy Positions
7. S.S Telescopic I.V Rod.
8. Detachable S.S Tray.

1,50,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

38

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.