

WS - 07

PLAIN BED

Size : 72" x 30" x 22"
 75" x 33" x 22"
 75" x 36" x 22"



M.S. Bow :- 8,000/-
 S.S. Bow :- 8,500/-
 M.S. Ply Bow :- 10,000/-
 S.S. Ply Bow :- 10,500/-

Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mm Dia. Optionally available with S.S. Side Rails for patient safety, Stainless Steel head and foot and bows, telescopic. Adjustable I.V.Rod, 100 mm dia. Revolving castors. Pretreated & epoxy powder Coated Supplied in KDC.

WS - 05

SEMI FOWLER BED

Size : 72"x36"x20"
 Size : 84"x35.5"x22"



M.S. Bow :- 9,500/- S.S. Bow :- 10,000/-
 M.S. Ply Bow :- 11,500/- S.S. Ply Bow :- 12,000/-

Made from uniformly perforated 18 Swg. CRCA Sheet mattress. Suitable for positions of Backrest with the help of a screw mechanism. Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mmDia. Optionally available with S.S. Side Rails for patient safety, telescopic adjustable I.V.Rod, 100mm dia. Revolving castors with wheels. Per treated & Epoxy powder coated, Supplied in KDC

WS - 04

FOWLER POSITION BED

Size : 81" x 35.5" x 22"

M.S. Bow :- 12,750/- S.S. Bow :- 13,250/-
 M.S. Ply Bow :- 15,250/- S.S. Ply Bow :- 15,750/-



Four section uniformly perforated 18 Swg. CRCA Sheet mattress. Suitable for positions of Backrest & Knee rest with the help of two screw mechanism. Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm x 30mm ERW Square tube Bows completely stainless steel 30mmDia. Optionally available with S.S. Side Rails for patient safety, telescopic adjustable I.V. Rod, 100 mm dia. Revolving castors with wheels, Per treated & epoxy powder coated, Supplied in KDC.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
 EQUIPMENTS PRIVATE LIMITED

Western Surgical

9

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.