


# WS - 07 E

## ECONOMY BED

Size : 75" x 32" x 22"



5,500/-

Made from 4mm Uniformly welded wiremesh with M.S. tubular head and foot end bows. Bed frame made from 50 x 25xE.R.W tube with Leggings made from 30x30mm Square ERW tube. Optionally available with S.S. Side Rails for patient safety, Stainless Steel head and foot and bows, telescopic. Adjustable I.V.Rod, 100 mm dia. Revolving castors. Pretreated & epoxy powder Coated. Supplied in KDC.



5,500/-

# WS - N 10

## ATTENDENT BED

Overall Size : 72" x 24" x 15" ht

Frame : 50mm x 20mm CRCA rectangular tube. Available with or without mattress. pre treated & epoxy powder coated. Mattress can be either built in or simple detachable.

# WS - 10

## DROP SIDE BABY BED

Overall Size : 54" x 30" x 22"

Bed with uniformly perforated sheet mattress 18 Swg, CRCA Sheet. With Safety side full length drop side. Pretreated & Epoxy Powder Coated. Supplied in KDC.



9,500/-

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

11

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.