



**WS - M5**

## **SEMI FOWLER POSITION BED MOTORIZED**

**Size : 2160 x 965 x 560 mm**

- Motorized with 2 Pcs actuators 1 Control box, 1 handset.
- Backrest & Footrest Motorized/ Manual
- Bed board in 16 SWG CRCA sheet
- Optional Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- Optionally available with battery backup

**56,500/-**

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

3

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.