



WS - M3



FIXED HELGHT ICU 5-FUNCTION SUPER DELUX BED MANUAL

Size: 2100 x 965 x 560mm

- Motorized with 2 Pcs actuators 1 Control box,
1 handset with battery backup
- Backrest & Footrest Motorized.
- Trendelenberg & Reverse Trendelenberg on gas Spring
- Bed board in High strength engineering plastic
- 4 underside occidental style bed side rails
- Wheels in high quality engineering plastic
- Head and foot end bows in high
strength engineering plastic
- 4 pcs uring hooks.
- 4 pole Saline attachment.

1,60,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

5

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.