



**WS - M8**

## ICU CLASSIC MANUAL DELUX

**Size : 2250 x 965 x 560 x 485-710mm**

- Manual adjustable 5 position bed with the help of high quality Screw mechanisms
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- 4 pole Saine attachment.
- 4 Saline Provisions with Aluminum inserts



52,500/-

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

6

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.