



WS - M10

FIXED HEIGHT ICU CLASSIC MANUAL

Size : 2150 x 965 x 555 mm

- Manual adjustable 4 position bed with the help of high quality Screw mechanisms
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- 4 pole Saine attachment.
- 4 Saline Provisions with Aluminum inserts



41,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

7

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.