

# Head Bow & Foot Bow Variants

We do not manufacture beds, we only deal in plastic components/spare parts



## MICHI - M1 HEAD/FOOT BOW

- Flagship product of the company.
- Robust, Ergonomic, Aesthetic, Smooth Finish
- Rod type fitment
- Center distance - **694mm** (+/- 2mm)

## MINI - M2 HEAD/FOOT BOW

- Flagship product of the company.
- Robust, Ergonomic, Aesthetic, Smooth Finish
- Rod type fitment
- Center distance - **850mm | 886mm** (+/- 2mm)



## MISAN - M3 HEAD/FOOT BOW

- Clip type & Hook type fitment
- Inbuilt buffer made from PC for shock absorption
- Side buffer has the provision to hold the IV rod
- Center distance - **Hook 695mm | Clip 750mm** (+/- 2mm)
- Design registration No : **342825-001**

## MIYON - M4 HEAD/FOOT BOW

- Clip type fitment
- Center distance - **910mm** (+/- 2mm)
- Design registration No : **328344-001**
- Patent registration No : **202021010270**



\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

16

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.