

# Head Bow & Foot Bow Variants

We do not manufacture beds, we only deal in plastic components/spare parts




## MIYON - M4 HEAD/FOOT BOW

- Clip type fitment for easy assembly
- Innovative concept of inbuilt CPR profile on the backside of the bow for **CPR emergencies**

## MIGO - M5 HEAD/FOOT BOW

- Rod type fitment
- Ideal for Deluxe, Semi-Fowler, Economical Beds
- Centre distance :  
**850mm | 869mm | 885mm | 900mm** (+/- 2mm)
- Design registration & patent awaited




## MIKU - M6 HEAD/FOOT BOW

- Robust, Ergonomic, Aesthetic, Smooth Finish Rod type fitment
- Centre distance : **694mm** (+/- 2mm)
- Design registration No : **357222-001**

## MIJU-M7 HEAD/FOOT BOW

- Elegant design, Eye-catching contours & profile, Good Finish, Strong body
- Centre Distance – **694 mm** (+/- 2mm)
- Design registration No : **372058-001**



## MIJUUNI - M8 HEAD/FOOT BOW

- Clip type fitment
- Ideal for Deluxe, Semi-Fowler, Economical Beds
- Centre distance : **910mm** (+/- 2mm)

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

17

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.