


Side Railing Variants

We do not manufacture beds, we only deal in plastic components/spare parts





MIROKU - MSR1 SIDE RAILINGS

- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator

MINANA - MSR2 SIDE RAILINGS

- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator
- Covers the full length of the bed
- Provision of embedded remote panel on side rail of MSR 2



MIHACHI - MSR3 SIDE RAILINGS


- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator
- Provision for urine bag holder
- Design registration No : 355150-001

MIJUSAN - MSR4 SIDE RAILINGS

- Duck down railing with bracket mechanism and provision of embedded remote panel
- Elegant design, Eye-catching contours and profile, Good Finish, Strong body
- Inbuilt membrane possible
- Push capacity of 80 kgs
- Provision of angle indicator
- Provision for urine bag holder
- Design registration No :

MSR4 Head Side: 372059-001

MSR4 Leg side: 372060-001



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

18

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.