


Mattress Platform

BED BOARD PLANK



BLOW MOULDED BED LINERS



- 4 section mattress platform for hospital beds
- X-ray compatible
- Easy to assemble
- Easy to clean / wash
- Lightweight but very sturdy
- Made out of Ultra High Density Poly Ethylene
- Modular system enables time saving during bed manufacturing
- Suitable for bed size 2270*942mm



Patient Report Holder

- A must need in every bed
- Made from High Density plastic
- Suitable for all types of beds
- Magnetic board on the front for writing patient's name
- Gaps available at the corners for easy cleaning
- Provision available for hanging on the head bow as well as wall mounted



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

19

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.