


TWO POINT IV BAG HOLDER MODEL NO.: SEP-IV01




DESCRIPTION

It is designed specifically for the overhead suspension of drip bags without the need to use valuable floor space. It is always available at a moment's notice, yet remains discreetly unobtrusive. As the carrier is track mounted, it can be moved to one side of the cubicle when not in use, and can be positioned at either side of the bed, or above it if required.

SALIENT FEATURES

- ◆ Strong aluminum track, white powder-coated finish
- ◆ Economical 2-point drip carriers with nylon runners
- ◆ Track suspended at cubicle track height
- ◆ Capacity 3 kg per hook
- ◆ Hooks made of stainless steel SS304
- ◆ Mobile track option
- ◆ Suspension system available in stainless steel S304 and MS powder coating.

TELESCOPIC IV BAG HOLDER MODEL NO.: SEP-IV02



DESCRIPTION

It is becoming increasingly important that the area around a hospital bed should be as uncluttered as possible. It is now accepted that the patient's quality of life can have a direct bearing on recovery time, and the perceived intrusiveness of medical equipment may be a contributory factor.

SALIENT FEATURES

- ◆ Heavy-Duty 5 Bag Holder
- ◆ Heavy duty self-locking wheeled trolleys
- ◆ Five point telescopic carrier trees
- ◆ Made of stainless steel SS304 and aluminum with white powder coating on it & hooks are of SS304
- ◆ Strong aluminum track with white powder - coated finish
- ◆ Track suspended at cubicle track height

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

24

© Copyright Reserved by : Western Surgical - Rajkot.