

WS- 1095A Folding Stretcher (canvas)

3,500/-

(213 L x 56 W x 15 H cms)
M.S. Pole width wise fold.



WS- 1095B Folding Stretcher 2 Fold

Unfolded Size : 221L x 53W x 15 Hcms
Folded Size : 110L x 18W x 11 Hcms
Net Wt : 5 kg
Gross Wt : 6 kg
Load Baring Capacity : 159 kg

5,500/-



WS- 1095C Folding Stretcher 4 Fold

Unfolded Size : 210L x 53W x 12 Hcms
Folded Size : 61L x 30W x 12 Hcms
Net Wt : 6 kg
Gross Wt : 7 kg
Load Baring Capacity : 159 kg

7,000/-



WS- 987A Scoop Stretcher (Telescopic/immobile)

A portable first aid stretcher, is commonly used for transporting bone fractured and seriously Injured patients. A clutch zip design at the middle of the model, the stretcher can be divided into left and right parts. Both parts can be drawn out when you need to put the patient on the operating table or the hospital bed quickly without moving the patient. The stretcher length can be directly transferred into the ambulance. The stretcher can fold or be separated into two parts, easy to carry or transport.

Technical Parameter :

Product size (LxWxH)	: 168x43x7cm(min)
	: 224x45x10cm(max)
Net Weight	: 9 kg
Gross Weight	: 11 kg
The bearing pressure	: < 159 kg
Dimension of the packing box	: 127x45x10cm



8,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

(27)

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.