



WS - 34

STRETCHER ON TROLLEY

Overall Size : 84" x 22" x 32"

Tubular steel structure with removable stretcher of M.S. Sheet Top. Optionally available with 4" or 7" dia. Wheels. Oxygen Cylinder Cage, Tray for Instruments, I.V. attachment with saline Rod. Pretreated & Epoxy Power Coated or completely Stainless Steel.

M.S. TOP :- 6,000/-
S.S. TOP :- 7,500/-

WS - 33A

SPECIAL STRETCHER ON TROLLEY

Overall Size : 84" x 22" x 32"

Tubular steel structure with removable stretcher of M.S. Sheet Top. available with 4" Whels. Oxygen Cylinder Coge, Tray for Instruments, I.V. attachment with saline Rod. Pretreated & Epoxy Power Coated



15,000/-

WS-34B

Stretcher on Trolley



WS-814

Patient Trolley (Hi-Low)

- ❖ Approx. size : 195L x 65W x 68-90H cms.
- ❖ Frame Work made of S. S. Tubes
- ❖ S. S. Side railing, I.V. rod.
- ❖ Adjustable Height by crank handle
- ❖ 15 cms dia castors.
- ❖ Full Stainless Steel.
- ❖ Also Available in complete M.S. frame.



25,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.