

# Patient Transfer Trolleys

9,500/-

USI- 991

## Patient transfer system

Patient mover from bed to bed

Applicable in large, medium and small hospitals, ICU, Operation Theatre, Emergency Room and Patient Ward, etc.

### Usage:

According to current requirement of modern nursing science, nursing involves a lot of physical work, which causes backache and lumber pain on parts of the nurses. While being used to move the patient, especially the patient who needs special care, between stretchers, patient beds, trolleys inspection beds, X-ray fluoroscopy tables and operation tables, EASYMOVE is easy to use.



### Technical Specifications:

Description	: Patient mover from bed to bed
Size	: 170 x 50 cm
Material	: Nylon, Plywood, Water-resistant.
Features	: Easy to be clean and sterilized. Maximum load for 15cm gap will be 200kg.



NURSES WORKS SIMPLER  
PATIENTS FEELS SAFE

## WS - 985 Wheelchair Stretcher with varied positions

35,000/-

### Technical parameter :

(Product size) (LxWxH): 190x59x85 cm  
 (Net weight): 45kg  
 (The bearing pressure): <159kg.  
 (Dimension of the packing box):  
 200x60x27cm  
 (Gross weight): 55kg



Automatic wheelchair type stretcher is mainly used to transport the sick, wooded outdoor. it's a necessary accessory on ambulance and could be folded as wheelchair in narrow stairs or elevators. This model is made of high quality aluminum alloy, with characteristics such as light weight, good appearances, reasonable configuration, convenience operation, antiseptis. Besides, it has soft foam matters block sidely. it's an easy operation stretcher.



## WS - 34Z

### TRAUMA CARE RECOVERY TROLLEY

Overall Size : 84" x 24" x 24-32"

- Height : Adjustable by foot operated Hydraulic pump
- X-Ray Permeable removable stretcher with built-in back rest
- Trendelenburg / Reverse trendelenburg position on gas spring
- Swing away Stainless steel side rails
- Oxygen Cylinder Holder
- Sliding X-Ray Cassette tray



65,000/-

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

31

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
 Ph. +91-281-2232871 / 2243159  
 Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.