



21,500/-

WS - 32

TRANSFER TROLLEY SYSTEM

"The exchange area technique has come to stay. This new concept prevents any possibility of cross infection. On the way to the operating suite the patient is brought to the exchange area on a ward trolley. The under carriage of the ward trolley and that of the operating room along together in the exchange area. With the trolley top the patient is transferred to the operating room undercarriage on easy rolling castors. The undercarriage is then separated and patient is taken to the operating room. The reverse process taken place after the operation when the patient is to be shifted to his ward.

Tubular steel structure with removable stretcher of M.S. Sheet Top. Available with 7" dia. Wheels. Oxygen Cylinder Cage, Tray for Instruments, I.V. Attachment with saline Rod. Pretreated & Epoxy Power Coated or completely Stainless Steel.

WS - 1008

Ambulance Stretcher

Light weight, strong structure, lovely design, good arts and crafts, easy antiseptis, small size, convenient to use.

Utilizes high precision mechanical transmission. Can be used automatically folded through the left and right hand control handle.

Made of high quality aluminum alloy tubes.

Dimensions :



33,000/-

Stretcher Trolley CT & MRI



45,000/-

A - High position : 1,900 x 560 x 900mm
 B - low position : 1,900 x 560 x 280mm
 Load weight : upto 120Kg
 Product weight : 40Kg
 Diameter of the castors : 125mm
 Applicable for various kind of hospitals, emergency centers, ambulances and battle fields

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

30

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.