

WS- 988 Spine Board

9,500/-



Spine Board is mainly used for immobilizing and safely transferring serious patients during rescuing. it's a first-aid transport equipments for medical department. it also could go into X-ray directly and float on water. This model is made is made of high strength engineering plastics, it can be enduring used, and has the advantage of antiaging, weight loading and easy stocking specialities.

The Features of our Spine Board are :

- CE certified Spine Board- in line with European Standard • Withstands up to 300 kgs.
- Available with 3 belts & groove to hold the board • Rigidity-Long life and light
- Light to carry • Floatation • Resistance to bump and corrosion • Does not absorb Liquids
- Immune to infiltration • Easy to clean (Water and soap are enough) • Translucent
- Homogeneity with X-Ray • Seamless moulding to eliminate ingress of fluid
- Avoids splintering and cracking
- Rigid Support During Movement of a Patient with Suspected Spinal or Limb Injuries., Ambulance Services and Lifeguards
- Better in Strength, User Friendly & Comply with Exacting Global Medical Standards

Specification of Spine Board

Sr. No.	Particulars	Standard	Test Method
1	Description	HMHDPE plastic moled stretcher for immobilization of patient.	Visual
2	Raw Material	Virgin HMHDPE granules form Lotrene Having grade TR-571 for Spine Board and Nylon for Belt.	Test Certificates
3	Colour	Orange	Visual
4	Total Article Wight (Kgs)	7.200 +/-2%	Electronic Weighing Balance
	Dimensional		
1	Article Length (mm)	1823 +/-5mm	Measuring Tape / Scale
2	Article Width (mm) Max.	448 +/-5mm	Measuring Tape / Scale
3	Article Thickness (mm)	61 +/-5mm	Measuring Tape / Scale
4	Load Capacity (Kgs.)	300.00	TM/QC/02

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

32

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.