

Delivery Table

Gynceologic Couch

Electric



1. Overall Size : 1900mm L x 620mm B x 700-950mm H .
2. Frame work made up of CRCA tabular pipe frame.
3. Three section cushioned mattress top with detachable leg section.
4. Trendlenburg / Reverse Trendlenburg Positions , Height Adjustment, Backsection adjustment through Actuator System.
5. Provided with a pair of Upholstery Lithonomy Positions
6. Detachable S.S Tray.

Remot :- 85,000/-
Hydrolic :- 65,000/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

37

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.