



WS - 45

INSTRUMENT TROLLEY S.S. TOP

Overall Size :	27" x 18" x 30"	30" x 20" x 30"
	4,500/-	6,000/-

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



WS - 47

DRESSING TROLLEY STEEL

Overall Size : 30" x 20" x 30"

Tubular structure with two stainless trays, having side rails with brackets for stainless bowl and bucket. Mounted on 100 mm dia. Revolving castors. Pretreated & Epoxy Power Coated or fully S.S.

9,000/-



INSTRUMENTS TROLLEY FULLY S.S. STEEL

Overall Size :	27" x 18" x 30"	30" x 20" x 30"
	7,500/-	9,500/-

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



WS - 245
RECT. STAINLESS
STEEL BED
SIDE TROLLEY

Overall Size
430 x 300 x 760mm ht.

6,500/-



WS - 246
SQ. STAINLESS
STEEL BED
SIDE TROLLEY

Overall Size
380 x 380 x 760mm ht.

6,800/-



21,000/-

Cathlab Trolley SS

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED