



INSTRUMENT TROLLEY S.S.

Overall Size : **27" x 18" x 30"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



INSTRUMENT TROLLEY S.S.

Overall Size : **30" x 20" x 30"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



INSTRUMENT TROLLEY S.S.

Overall Size : **48" x 24"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



3 SELF STEEL TROLLEY

Overall Size : **20" x 30"**

A-Steel Top
B-Fully Steel



3 SELF STEEL TROLLEY

Overall Size : **18" x 15"**

A-Steel Top
B-Fully Steel

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

40

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.