

Any Requirements Sizes Are Also Available We Manufacturing

WS-38 OVER BED TABLE PUSH TYPE

44" x 17.5" x 38" ht

Tubular Steel Structure with laminated Tray on 50mm dia, revolving castors pre-treated & powder coated



4,100/-

WS - 36

7,200/-

ADJUSTABLE BED SIDE TABLE

Overall Size : 24" x 18" x32" - 50"

Height adjustable, tubular steel structure with laminated top, Adjustable on S.S. Telescopic tube. Mounted on 50mm dia. Revolving castors. Pretreated & Epoxy Powder Coated. Also Available in Fully Stainless Steel. Optional available fully s.s.



ADJUSTABLE BED SIDE TABLE GEAR TYPE

Overall Size : 30" x 16" x32" -50"

Height adjustable with help of gear mechanism provided on the side of the table. Tubular steel structure with laminated top with adjustable tilt. Mounted on 50mm dia. Revolving Castors. Pre treated & Epoxy Powder Coated. Also Available in Fully Stainless Steel.

7,200/-

WS - 36 A



Available in Attractive Colours



WS 601

OVERBED TABLE FOLDING WITH HYDRAULIC FUNCTION



6,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

45

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.