

Scrub Station & Instrument Washing



Dimension's As Per Your Size



Dimension's As Per Your Size



4 Feet

Dimension's As Per Your Size



65,000/-



Dimension's As Per Your Size

WS - X13

SCRUB SINK

Size : 20" d X 39" h X vertical depth 18"

Scrub sink made in complete 18 swg 304 grend SS Sheet. Tailor made to fit any length, minimum 30".

Drainage options either on left or right Available with or without bottom cover Can be provided with foot pedal taps automatic sensor taps. Soap Dispensers

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

46

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.