

WS - 13

SOILED
LENIN BIN

Overall Size : 20" dia X 36" ht

Tubular Structure mounted on 100mm dia revolving castors. Provided with detachable canvas back. Pre-treated epoxy powder coated / Complete Stainless Steel.



2,500/-

WS - 15 G

OPERATION TABLE FOR
ANIMAL SURGERY

HEIGHT ADJUSTABLE

Size : 48" x 26" x 30" ht.

Epoxy Powder Coated. Completely Stainless steel



11,000/-

KICK BUCKET SMALL



2,100/-

WS - 54

WS - 57

BOWL STAND

Overall Size : 33" X 14" dia with
Tubular Structure mounted on 50mm dia. revolving castor with two stainless steel bowl. Pretreated & Epoxy Powder Coated / completely SS.



2,700/-

WASH BASIN SINGLE

Overall Size : 35" X 14" dia with
Tubular Structure mounted on 50mm dia. revolving castor with a stainless steel bowl. Pretreated & Epoxy Powder Coated / completely SS.



2,600/-



3,600/-

**Wash Basin
MS & SS**

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.