

Any Requirements Sizes Are Also Available We Manufacturing



3,600/-
Fully Steel



1,800/-
S. S. Top



1,500/-
S. S. Top



1,500/-

WS - 39

REVOLVING STOOL

Overall Size : 11.5" dia x 21" - 29"

Tubular Steel Structure with Superimposed Stainless Steel Top. Height adjustable from 425 to 675 mm with the help of smooth screw mechanism. Pretreated & Epoxy Power Coated. Oponationally available fully Stainless Steel.



4,100/-
Hydrolic with Wheel



3,000/-
Hydrolic Fix



S. S. Top
1,500/-



M.S. Top
1,200/-



Fully Steel
2,700/-



WS 393 : REVOLVING STOOL WITH CUSHION.
Size : 370 dia 545-745mm ht. / 14" dia 17"-29"
Available in MS Powder coated / S. S 304 grade.



WS 393 : REVOLVING STOOL WITH CUSHION.
Size : 370 dia 545-745mm ht. / 14" dia 17"-29"
Available in MS Powder coated / S. S 304 grade.

WS - 68

BED SIDE STOOL (Overall Size : 12" x 12" x 8")

Mode from CRCA steel structure with M.S. Pretreated Epoxy Power Coated. Oponationally available with laminated top, Stainless Steel top or fully Stainless Steel.



3,200/-

Food Trolley on Bed
Size : 18" x 24" x 15" H.
Complete Steel & Powder Coted

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

47

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.