



15,000/-

**WS - 66**

**SURGICAL INSTRUMENT CABINET**

**Size : 78" X 19.5" d X 36" W**

Compact Cabinet for instruments provided with 4 shelves to form 5 compartments. Front doors in glass 5mm thick.



18,000/-

**WS - 67**

**INSTRUMENT CABINET**

**Size : 63" X 17.5" d X 25" W**


Frame : CRCA rectangular tube, CRCA sheet plants  
Five Shelves, door & two sides of 5mm thick fine glass. Pre-treated & epoxy powder coated.

**Surgical Drawer**

**Surgical Cupboard**

**Multi Locker - 4**

**Multi Locker - 12**



15,000/-



22,000/-



8,500/-  
10,500/-



**Multi Locker - 12**

**WS 673 : CABINET 6 COM.**

**WS 674 : CABINET 15 COM.**



Size : 1970 x 500 x 920mm ht.  
78" x 19" x 36"  
0.8mm thick CRCA sheet  
Completely Powder Coated.



Size : 1970 x 500 x 920mm ht.  
78" x 19" x 36"  
0.8mm thick CRCA sheet  
Completely Powder Coated.

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

50

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.