



WS-

4,500/-



WS-

5,000/-



4,500/-



6,500/-

Surgeon Stool

Tubular Steel Structure with Superimposed Stainless Steel Top. Height adjustable with the help of Heavy Duty gas spring. Back made up of uniformly perforated CRCA. Sheet. Pretreated & Epoxy Power Coated.



WS- 1055

Surgeon Chair

Overall Approx. Size : Cushioned seat and back
Pneumatic height adjustable
Height adjustable and swinging back rest
5 Legs chrome plated base with foot ring
Base fitted with 50 mm dia castors

WS - 59

Overall Size :
20"x18"x18"x20"x12"x8"
Tubular Structure with M.S.
Top with Aluminum
Treads. Pretreated
and Epoxy Powder
Coated.



Single Step
Stool M.S.



Double Step
Stool M.S.



Single Step
Stool Fully Steel



Double Step
Stool Fully
Steel

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

52

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.