

Economy Laparoscopy / Monitor Trolley



Special Features

1. Size - 4.75 Ft. H x L + W - 2.25 Feet
2. Attractive Color Combination
3. Self Size - 18" x 20" tray type
4. Self Height Adjustable No. 2,3,4
5. One Drawer with Lock
6. 4" Castor with Lock systems
7. Two M. S. Handle for Movement
8. Co2 Cylinder stand Back Side with 2 clamp
9. 6 Plug 15 x 5 Amp Switch with 10 Ft. Cable
10. Space for UPS / Servo Stabilizer behind



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

56

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.