

## WS- 939

### Emergency Trolley (ABS)

ABS body, metal ball-bearing rail. A-Aluminium column B- Plastic steel columns. IV Pole, needle disposal holder (optional); Dust basket, utility container, file bag, defibrillator shelf, power outlet & hooks, Oxygen tank holder (Optional), C.P.R. Board (Optional), Sliding side shelf, stainless steel guard rail, concave top covered transparent soft plastic glass, Centralized lock, 2 small & 3 big drawers repositioned easily with partitions can be well organized. ABS basket & drawer frame can be repositioned easily. Luxurious wise casters, two with brakes.

Approximate Size : 850 x 520 x 920 mm

ON DEMAND



## WS- 940

### Medicine Trolley (Cassette Type)

- \* Designed to supply medicine when nurse make around of visit to patient
- \* Two sides unit cassette and one side plate cassette
- \* Easy to discharge
- \* Installed multi-bin container at rear side to make easy and multi possible strong
- \* Multi structure can be possible used by kinds of accessories
- \* Ergonomics designed push-handle and efficient structure
- \* Total lock key system

Size : 780(L) x 490 (D) x 920(H)



ON DEMAND

## WS- 40D

### Case History Trolley

Approx Size :  
680 x 400 x 850 mm

Made of ABS Material

Double row, 40 layers  
with a drawer

3" Castor

ON DEMAND



## WS- 20D

### Case History Trolley

Approx Size :  
336 x 400 x 850 mm

Made of ABS Material

Single row, 20 Layers  
with a drawer

3" Castor



\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.