

POWER COLUMN

Power Column samples the organization of medical service equipment in one aesthetically appealing vertical headwall. It is a modular designed Product to provide a clean, sanitary and ordered environment. Furthermore it can enhance the quality in care and medical environment. Due to its systematized design, a network provided for medical care service is applicable to the use of secondary instruments and assemblies. The main structure consists of fluorescent panel, aluminum component and side track. Also includes power lines and medical gas piping. Side tracks are installed at four at angles. and the length is equal to the column. All the assemblies required by medical care personnel can be easily and securely placed on the side track of this column. Applications in mid to acute level environment



Features

- ❖ Designed for ease of maintenance and infection control.
- ❖ Built as a four-sided unit that offers 360° access to services and equipment management for critical care areas.
- ❖ Can be customized and upgraded to meet your facility's needs.
- ❖ Offered in a variety of widths and heights, tailored to meet your space needs.
- ❖ Wide selection of laminate colors are available.
- ❖ Multiple rail system allows for superior accessory management.
- ❖ Removable laminate panels allow for simplified maintenance.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

72

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.