


ONE STOP SHOWROOM FOR HOSPITAL & HEALTHCARE PRODUCTS.



HOSPITAL FURNITURE

CATALOGUE - 2024



Western Surgical

®

WESTERN SURGICAL EQUIPMENTS PRIVATE LIMITED

AN ISO 9001 : 2008 & ISO 13485 : 2015 AC 2012, CE, FDA & GMP CERTIFIED COMPANY

www.westernsurgical.in Email : westernsurgical@gmail.com

Our Valuable Clients

Sr. No.	Dr. Name & Hospital Name	City Name
1	Dr. Viral Vasavada (general Surgeon) M. S. Samarpan Hospital	RAJKOT
2	Dr. Sanjay Pandya M.D., D.N.B., Nephro	RAJKOT
3	Dr. Raju Sagar M.S. Laparoscopy Surgeon	RAJKOT
4	Dr. I.P. Swami M.D. Physician	RAJKOT
5	Dr. Pradip Shah Pediatric	RAJKOT
6	Dr. Rajesh Gandhi Orthopedic	RAJKOT
7	Dr. Rakesh Chokshi Laproscopy Surgeon	RAJKOT
8	Dr. Sanjay Patoliya M.S. Bariatric Surgeon	RAJKOT
9	Dr. Shrenik Shah M.S. ENT	RAJKOT
10	Dr. Himanshu Thakkar M.S. ENT	RAJKOT
11	Dr. Samir Thakkar Pediatric	RAJKOT
12	Dr. Vinay Majmudar M.D. Physician	RAJKOT
13	Dr. Hitesh Dhru M.S. MCH. Plastic	RAJKOT
14	Dr. Girish Amlani M.S. MCH. Plastic	RAJKOT
15	Dr. Bhaumik Bhayani M.S. MCH. Plastic	RAJKOT
16	Dr. Ashish Jasani M.S. MCH. Urologist	RAJKOT
17	Dr. Bhavesh Sachde Orthopedic Join Replacement Surgeon	RAJKOT
18	Guru Govindsinh Hospital General Hospital	JAMNAGAR
19	Civil Hospital - RAJKOT, AMRELI, SURAT, AHMEDABAD, MORBI, PORBANDAR, ETC.	
20	Aashirvad Multispicillity Hospital	RAJKOT
21	Bhavsinhji General Hospital	PORBANDAR
22	Sir-t Hospital General Hospital	BHAVNAGAR
23	Swami Nirdoshanand Hospital - Timbi Multispeciality Hospital	TIMBI
24	Samarpan General Hospital	JAMNAGAR
25	C. U. Shah Medical College	SURENDRANAGAR
26	Christ Hospital General Hospital	RAJKOT
27	Sadbhavna Hospital	MAHUVA
28	City Dental Hospital	RAJKOT
29	Giriraj Multispeciality Hospital	RAJKOT
30	H. J. Doshi Sarvajanic Hospital General Hospital	RAJKOT

For More Valuable Customer on our web site www.westernsurgical.in

INDEX

ICU Excel Motorized.....	1
Fixed Height ICU Classic Motorized.....	2
Semi Fowler Position Bed Motorized.....	3
Fowler Position Bed Classic Motorized.....	4
Fixed Height ICU 5- Function Super Delux Bed Manual.....	5
ICU Classic Manual Delux.....	6
Fixed Height ICU Classic Manual.....	7
Plain Bed Excel / Semi Fowler Bed Excel	
Fowler Position Bed Excel.....	8,10
Economy Bed / Attendent Bed / Drop Side Baby Bed.....	11
Bed Illustration with Hospital Furniture Parts.....	12,18
Mattress Platform & Patient Report Holder.....	19
Mattress for Hospital Beds	20
Bed Side Riling.....	21
Twin Bed For Medical Hostel.....	21
IV Stand	22,24
Bed Side Locker.....	25
Oxygen Cylinder Trolley.....	26
Stretcher.....	27,31
Spine Board.....	32
Waste Carrying Trolley.....	33
Bio Medical Waste Bins.....	33
Moping Trolley.....	33
BMW Box Type Trolley.....	33
CSSD Air Tight Trolley.....	33
Examination Table & Gynec Table & Physician Table.....	34
Operation Table Screw Type & Delivery Table.....	35
Telescopic Delivery Table.....	36
Attendance Bed.....	36
Delivery Table.....	37,38
Instrument Trolley.....	39,40
ECG Trolley & Derma Trolley & Defibrillator Trolley.....	42,43
Laundry Trolley & Retina Machine Trolley.....	43
Fresh Soild Lenin Trolley.....	43
Mayos Hight Adjustable Trolley & Table	44
Laser Machine Trolley.....	44
V. H. Trolley.....	44
Over Bed Table Push Type & Adjustable Bed Side Table.....	45
Over Bed Table Folding with Hydraulic.....	45
Scrub Station & Instrument Washing.....	46
Revolving Stool.....	47
Food Trolley on Bed.....	47
Soiled Lenin Bin.....	48
Operation Table Animal Surgery.....	48
Kick Bucket Small & Bowl Stand & Wash Basin Single.....	48
Waiting Chair.....	49

INDEX

Surgical Instrument Cabinet.....	50
Bed Side Screen.....	51
Surgeon Stool.....	52
Single & Double Step Stool.....	52
Emergency Crash Cart Trolley.....	53
Emergency I.C.U. Defibrillator Trolley.....	54
Emergency Drub Trolley.....	54
Trauma Care Medicine Trolley.....	54
Endo & Laparoscopy Trolley.....	55
Economy Laparoscopy & Monitor Trolley.....	56
Perforated E.T.O. Tray & Net Type E.T.O. Tray.....	57
Cidex Tray for Lapro & Urology.....	57
Burn Care & Catheter Hook Cardiac OT Lab.....	58
Variable Height Tilt Table.....	58
Drug Trolley Twin.....	59
X-Ray Apron Stand & Patient Shifting Rollar.....	59
Syringe Pump Stand.....	59
Emergency Trolley.....	60
Medicine Trolley (Cassette Type).....	60
Case History Trolley.....	60
Blood Collection Pathology Chair.....	61
Blood Transfusion Chair Delux.....	61
Attendent Bed & Chair.....	62
Sofa Chair Twin Seater.....	62
Sofa Chair Single Seater.....	62
Quadriceps Exercise Table.....	63
Shoe Rack Steel.....	65
Shoe Rack Steel.....	65
Medical Rack Surgical Medicine File.....	65
Medical Student Chair.....	66
Medical Equipments & Instrument Safety Locker.....	67
Doctor's Chair.....	68
Doctor Office Table.....	68
Medical Computer Table.....	68
Hospital Visitor Chair.....	69
Customized Plastic Moulded Chair Parts.....	70
Curtain Tracks.....	71
Power Column.....	72
Surgical / Anaesthesia Pendant.....	73
Single Arm Surgical / Anaesthesia Pendant.....	74
Wall / Ceiling Mounted I.C.U. Pendant.....	74
Under Bed Crib.....	75
Crib On Stand.....	75
Accessories Hospital Furniture.....	76



WS - M1



Size 2150 x 965 x 510 - 730mm ht.

ICU EXCEL MOTORIZED

- Completely motorized 4 Pcs actuators 1 Control box, 1 Pc handset with battery backup
- Bed board in High strength engineering plastic
- Wheels in High quality engineering plastic With central Locking
- Head and foot end bows in high strength engineering plastic

- 4 underside occidental style bed side rails
- 4pcs urine hooks.
- 4 pole Saline attachment.
- Optionally available with urgent release (CPR) function
- Also Available with ACP if required

1,75,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

①

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M4



**FIXED HELGHT ICU
CLASSIC MOTORIZED**

Size: 2100 x 965 x 460mm

- Motorized with 2 Pcs actuators 1 Control box, 1 handset with battery backup
- Backrest & Footrest Motorized.
- Trendelenberg & Reverse Trendelenberg on gas Spring
- Wheels in high quality engineering plastic
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs uring hooks.
- 4 pole Saline attachment.

1,35,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M5

SEMI FOWLER POSITION BED MOTORIZED

Size : 2160 x 965 x 560 mm

- Motorized with 2 Pcs actuators 1 Control box, 1 handset.
- Backrest & Footrest Motorized/ Manual
- Bed board in 16 SWG CRCA sheet
- Optional Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- Optionally available with battery backup

56,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

3

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M5

FOWLER POSITION BED CLASSIC MOTORIZED

Size : 2160 x 965 x 560 mm

- Motorized with 2 Pcs actuators 1 Control box, 1 handset.
- Backrest & Footrest Motorized/ Manual
- Bed board in 16 SWG CRCA sheet
- Optional Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- Optionally available with battery backup

56,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

4

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M3



FIXED HELGHT ICU 5-FUNCTION SUPER DELUX BED MANUAL

Size: 2100 x 965 x 560mm

- ▶ Motorized with 2 Pcs actuators 1 Control box,
1 handset with battery backup
- ▶ Backrest & Footrest Motorized.
- ▶ Trendelenberg & Reverse Trendelenberg on gas Spring
- ▶ Bed board in High strength engineering plastic
- ▶ 4 underside occidental style bed side rails
- ▶ Wheels in high quality engineering plastic
- ▶ Head and foot end bows in high
strength engineering plastic
- ▶ 4 pcs uring hooks.
- ▶ 4 pole Saline attachment.

1,60,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

5

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M8

ICU CLASSIC MANUAL DELUX

Size : 2250 x 965 x 560 x 485-710mm

- Manual adjustable 5 position bed with the help of high quality Screw mechanisms
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- 4 pole Saine attachment.
- 4 Saline Provisions with Aluminum inserts



52,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

6

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - M10

**FIXED HEIGHT ICU
CLASSIC MANUAL**

Size : 2150 x 965 x 555 mm

- Manual adjustable 4 position bed with the help of high quality Screw mechanisms
- Bed board in 16 SWG CRCA sheet
- Imported collapsible bed side rails
- Head and foot end bows in high strength engineering plastic
- 4 pcs urine hooks.
- 4 pole Saine attachment.
- 4 Saline Provisions with Aluminum inserts



41,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

7

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



WS - 701

PLAIN BED EXCEL

Size : 1900 x 900 x 560mm
1900 x 825 x 560mm

17,000/-



WS - 501

**SEMI-FOWLER
BED EXCEL**

Size : 2055 x 900 x 560mm

15,000/-



WS - 401

**FOWLER POSITION
BED EXCEL**

Size : 2055 x 900 x 560mm

17,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

8

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

WS - 07**PLAIN BED**

Size : 72" x 30" x 22"
 75" x 33" x 22"
 75" x 36" x 22"



M.S. Bow :- 8,000/-
 S.S. Bow :- 8,500/-
 M.S. Ply Bow :- 10,000/-
 S.S. Ply Bow :- 10,500/-

Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mm Dia. Optionally available with S.S. Side Rails for patient safety, Stainless Steel head and foot and bows, telescopic. Adjustable I.V.Rod, 100 mm dia. Revolving castors. Pretreated & epoxy powder Coated Supplied in KDC.

WS - 05**SEMI FOWLER BED**

Size : 72"x36"x20"
 Size : 84"x35.5"x22"



M.S. Bow :- 9,500/- S.S. Bow :- 10,000/-
 M.S. Ply Bow :- 11,500/- S.S. Ply Bow :- 12,000/-

Made from uniformly perforated 18 Swg. CRCA Sheet mattress. Suitable for positions of Backrest with the help of a screw mechanism. Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mmDia. Optionally available with S.S. Side Rails for patient safety, telescopic adjustable I.V.Rod, 100mm dia. Revolving castors with wheels. Per treated & Epoxy powder coated, Supplied in KDC

WS - 04**FOWLER POSITION BED**

Size : 81" x 35.5" x 22"



M.S. Bow :- 12,750/- S.S. Bow :- 13,250/-
 M.S. Ply Bow :- 15,250/- S.S. Ply Bow :- 15,750/-

Four section uniformly perforated 18 Swg. CRCA Sheet mattress. Suitable for positions of Backrest & Knee rest with the help of two screw mechanism. Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm x 30mm ERW Square tube Bows completely stainless steel 30mmDia. Optionally available with S.S. Side Rails for patient safety, telescopic adjustable I.V. Rod, 100 mm dia. Revolving castors with wheels, Per treated & epoxy powder coated, Supplied in KDC.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
 EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54



WS - 07

PLAIN BED

M.S. Ply Bow :- 10,000/-
S.S. Ply Bow :- 10,500/-

Size : 72" x 30" x 22" 75" x 33" x 22" 75" x 36" x 22"

Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mm Dia. Optionally available with S.S. Side Rails for patient safety, Stainless Steel head and foot and bows, telescopic. Adjustable I.V.Rod, 100 mm dia. Revolving castors. Pretreated & epoxy powder Coated Supplied in KDC.



WS - 05

SEMI FOWLER BED

**Size :
81" x 35.5" x 22"**

M.S. Ply Bow :- 11,500/-
S.S. Ply Bow :- 12,000/-



Made from uniformly perforated 18 Swg. CRCA Sheet mattress. Suitable for positions of Backrest with the help of a screw mechanism. Made from 18 Swg. Uniformly perforated sheet mattress Bed frame made in 60 x 30mm Rect. ERW tube. Leggings made in 30mm X 30mm ERW Square tube Bows completely stainless steel 30mmDia. Optionally available with S.S. Side Rails for patient safety, telescopic adjustable I.V.Rod, 100mm dia. Revolving castors with wheels. Per treated & Epoxy powder coated, Supplied in KDC

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54


10

© Copyright Reserved by : Western Surgical - Rajkot.

WS - 07 E

ECONOMY BED

Size : 75" x 32" x 22"



5,500/-

Made from 4mm Uniformly welded wiremesh with M.S. tubular head and foot end bows. Bed frame made from 50 x 25x E.R.W tube with Leggings made from 30x30mm Square ERW tube. Optionally available with S.S. Side Rails for patient safety, Stainless Steel head and foot and bows, telescopic. Adjustable I.V.Rod, 100 mm dia. Revolving castors. Pretreated & epoxy powder Coated. Supplied in KDC.

WS - N 10

ATTENDENT BED

Overall Size : 72" x 24" x 15" ht

Frame : 50mm x 20mm CRCA rectangular tube. Available with or without mattress. pre treated & epoxy powder coated. Mattress can be either built in or simple detachable.



5,500/-

WS - 10

DROP SIDE BABY BED

Overall Size : 54" x 30" x 22"

Bed with uniformly perforated sheet mattress 18 Swg, CRCA Sheet. With Safety side full length drop side. Pretreated & Epoxy Powder Coated. Supplied in KDC.



9,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159


Mo. +91 94280 03750/52/54

11

© Copyright Reserved by : Western Surgical - Rajkot.


HEAD BOW/FOOT BOW

These panels are the flagship product of the company. The panels are robust, ergonomic, aesthetic and come with a finish that complies with customer needs. For an unparalleled convenience of use, with help of connecting clamps & rods, they are extremely easy to assemble and remove. There are a variety of sticker color options available for these panels to enhance the aesthetics and sync with the interiors of the hospitals.



SIDE RAILING

Side railings complements the Head & Foot bow aesthetically. With a curved finish of our side rail, it enhances the look of the whole bed. There is a provision of angle indicator on the side rail which helps knowing the exact position of the bed while using motorized beds. Two side rails on one side covers the full length of the bed, thus not leaving any trap gaps in the bed.



BRACKET MECHANISM

Our side railing gives the feature of duck down railing with just a push of a button with the help of gas spring in the bracket. The push capacity is 80 kgs. We have an in-house manufacturing facility of these brackets with separate assembly line for the fitments.

BED BOARD PLANKS/BED LINERS

We have innovated the conventional mattress platform in plastic with two varieties 10-piece plastic bed board and 4-piece bed board to improve the bearing capacity of the bed board. The strength is much higher than the metal one as it shares the load equally. The 10 pc bed board plank is easy to attach with screws. The 4 pc bed liner is also easy to assemble on the metal bed frame. It is MRI & X-Ray compatible too. Both models are extremely time saving and also reduces the powder coating activity with respect to cost and labor.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *


Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

12

© Copyright Reserved by : Western Surgical - Rajkot.



FITMENT TYPE

We have developed a variety of fitments for mounting bow to the bed frame easily. These are developed keeping in mind the functional requirements & are made of a combination of materials that makes them strong, sturdy, presentable & user friendly. Plastic rods fitment specially used for M1 & M2 in place of metal ones which are non-corrosive & user-friendly. We also have developed Hook & Clip metal fitments, specially used for M3 & M4.

BUFFER WHEELS

These plastic buffers work as great substitutes for rubber / PVC buffers. It provides safety from wall surface score marks on the head bow. They have good shock absorption property. Available in customized colors & sizes. Outer dia - 5 Inch & inner dia 10.2, 14.5, 26.4 and 32.6 mm options. There is a special type of buffer used in M3. It has a provision to hold the IV Rod pole.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical


* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

13

© Copyright Reserved by : Western Surgical - Rajkot.

Bed illustrations with hospital furniture parts

We do not manufacture beds, we only deal in plastic components/spare parts




MICHU - M1 BOW BED

- Curved head /foot bow - M1 with rod type fitment
- MSR - 2 side railings with hydraulic bracket mechanism
- Side buffer wheel for safety from wall surface scar marks on the bed

MINI - M2 BOW BED

- Straight Head /Foot Bow - M2 With Rod Type Fitment
- MSR - 1 side railings with hydraulic bracket mechanism
- Side buffer wheel for safety from wall surface scar marks on the bed




MISAN - M3 BOW BED

- Straight head/foot bow - M3 with inbuilt buffer for shock absorption
- Clip type & hook type fitment
- MSR - 2 Side Railings With Hydraulic Bracket Mechanism

MIYON - M4 BOW BED

- Straight Head /Foot Bow - M4 With CPR Board Provision
- Clip type fitment
- MSR - 1 side railings with hydraulic bracket mechanism



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.




WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Bed illustrations with hospital furniture parts

We do not manufacture beds, we only deal in plastic components/spare parts





MIGO - M5 BOW BED

- Grilled head/foot bow - M5 with rod type fitment

MIKU - M6 BOW BED

- Curved head/foot bow - M6 with rod type fitment
- MSR-3 side railings with hydraulic bracket mechanism




MIJU- M7 BOW BED

- Curved head/foot bow - M7 with rod type fitment
- MSR-4 Side Railing with bracket mechanism and embedded panel

MIJUUNI - M8 BOW BED

- Clip type fitment
- Ideal for Deluxe, Semi-Fowler, Economical Beds



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Head Bow & Foot Bow Variants

We do not manufacture beds, we only deal in plastic components/spare parts




MICHI - M1 HEAD/FOOT BOW

- Flagship product of the company.
- Robust, Ergonomic, Aesthetic, Smooth Finish
- Rod type fitment
- Center distance - **694mm** (+/- 2mm)

MINI - M2 HEAD/FOOT BOW

- Flagship product of the company.
- Robust, Ergonomic, Aesthetic, Smooth Finish
- Rod type fitment
- Center distance - **850mm | 886mm** (+/- 2mm)



MISAN - M3 HEAD/FOOT BOW

- Clip type & Hook type fitment
- Inbuilt buffer made from PC for shock absorption
- Side buffer has the provision to hold the IV rod
- Center distance - **Hook 695mm | Clip 750mm** (+/- 2mm)
- Design registration No : **342825-001**

MIYON - M4 HEAD/FOOT BOW

- Clip type fitment
- Center distance - **910mm** (+/- 2mm)
- Design registration No : **328344-001**
- Patent registration No : **202021010270**



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Head Bow & Foot Bow Variants

We do not manufacture beds, we only deal in plastic components/spare parts




MIYON - M4 HEAD/FOOT BOW

- Clip type fitment for easy assembly
- Innovative concept of inbuilt CPR profile on the backside of the bow for **CPR emergencies**

MIGO - M5 HEAD/FOOT BOW

- Rod type fitment
- Ideal for Deluxe, Semi-Fowler, Economical Beds
- Centre distance : **850mm | 869mm | 885mm | 900mm** (+/- 2mm)
- Design registration & patent awaited




MIKU - M6 HEAD/FOOT BOW

- Robust, Ergonomic, Aesthetic, Smooth Finish Rod type fitment
- Centre distance : **694mm** (+/- 2mm)
- Design registration No : **357222-001**

MIJU-M7 HEAD/FOOT BOW

- Elegant design, Eye-catching contours & profile, Good Finish, Strong body
- Centre Distance – **694 mm** (+/- 2mm)
- Design registration No : **372058-001**



MIJUUNI - M8 HEAD/FOOT BOW

- Clip type fitment
- Ideal for Deluxe, Semi-Fowler, Economical Beds
- Centre distance : **910mm** (+/- 2mm)

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com


Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

Side Railing Variants

We do not manufacture beds, we only deal in plastic components/spare parts



MIROKU - MSR1 SIDE RAILINGS

- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator



MINANA - MSR2 SIDE RAILINGS

- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator
- Covers the full length of the bed
- Provision of embedded remote panel on side rail of MSR 2



MIHACHI - MSR3 SIDE RAILINGS

- Duck down railing with bracket mechanism
- Push capacity of 80 kgs
- Provision of angle indicator
- Provision for urine bag holder
- Design registration No : 355150-001



MIJUSAN - MSR4 SIDE RAILINGS

- Duck down railing with bracket mechanism and provision of embedded remote panel
- Elegant design, Eye-catching contours and profile, Good Finish, Strong body
- Inbuilt membrane possible
- Push capacity of 80 kgs
- Provision of angle indicator
- Provision for urine bag holder
- Design registration No :

MSR4 Head Side: 372059-001

MSR4 Leg side: 372060-001

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

18


© Copyright Reserved by : Western Surgical - Rajkot.

Mattress Platform

BED BOARD PLANK



BLOW MOULDED BED LINERS



- 4 section mattress platform for hospital beds
- X-ray compatible
- Easy to assemble
- Easy to clean / wash
- Lightweight but very sturdy
- Made out of Ultra High Density Poly Ethylene
- Modular system enables time saving during bed manufacturing
- Suitable for bed size 2270*942mm



Patient Report Holder

- A must need in every bed
- Made from High Density plastic
- Suitable for all types of beds
- Magnetic board on the front for writing patient's name
- Gaps available at the corners for easy cleaning
- Provision available for hanging on the head bow as well as wall mounted



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

19

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

WS - 5001

Mattress for
Hospital Beds




A-Single Mattress for plain beds.



**B-Two section mattress for
Semi-Fowler beds.**

**C-Four section
mattress for Fowler
I.C.U. Beds.**



- 1) Mattress for Attendance bed Size
- 2) Mattress for Exam Table Plain Size
- 3) Mattress for Exam Table Gynec Cut Size
- 4) Mattress for Labour Table Size 2.5x6
- 5) Mattress for Labour Table 3 Section Size
- 6) Mattress for Pediatric Bed Size
- 7) Mattress for Hydrolic ot Table Size
- 8) Mattress for Reddle Size
- 9) Mattress for Strecher Trolley Size

3,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(20)


© Copyright Reserved by : Western Surgical - Rajkot.

**WS R2 :
SLIDE DOWN
BED SIDE RAIL**



**Size : 1370 x 460mm ht.
54" x 18" ht.**

**Bed Side Reling
S.S. Allu.
Folding Imp.**



6,500/- Pair + GST

**Steel Collapsible Side Rails
(Set of 2) Approx Wt : 5.500 Kgs per Set**



10,500/-

Twin Bed for Medical Hostel

WS - 5004

Mosquito Net Poles

Four vertical M.S. tubes of 19mm,
Two longitudinal & Two transverse
tubes of 15.8mm.

Knocked down condition.

Pre-treated and Epoxy Powder Coated.



3,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(21)

© Copyright Reserved by : Western Surgical - Rajkot.



**CASTING
BASE
Delux**

1100/-



**CASTING
BASE**

900/-



**FIBRE
BASE
5 Wheel**

3,000/-



**STEEL
5 Wheel**

3,000/-



**M.S.
5 Wheel**

3,600/-



5 Wheel



850/-

WS - 28

SALINE STAND

Tubular structure with diepressed five pronged legs chrome plated rod. Pretreated and Epoxy Power Coated. Optionally, available with or without 50 mm dia. Revolving castors. Also available Saline Stand completely S.S. with fibre base.



550/-

**Ceiling
I. V. Stand
5 Hook**



3,600/-

On Bed I. V. Stand Steel

On Bed I. V. Stand MS

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.




**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

TELESCOPIC IV BAG HOLDER

(MODEL NO. : SEP-IV03)



DESCRIPTION

We offer you a complete line of track mounted IV bags holders with 5-bag pendants. The holders suspend from an Eze-Lock IV Carriage, and are made from stainless steel with folding arms for easy sterilization and storage.

SALIENT FEATURES

- ◆ Heavy-Duty 5 Bag Holder
- ◆ Heavy duty self-locking wheeled trolleys
- ◆ Five point telescopic carrier trees
- ◆ Made of stainless steel SS304
- ◆ Hooks are of SS 304.
- ◆ Strong aluminum track with white Powder-coated finish
- ◆ Track suspended at cubicle track height
- ◆ IV hanger closed length: 24 inches and openable length:36 inches with trolley

HOSPITAL CURTAINS



DESCRIPTION

We offer an economical range of hospital cubical curtains. All varieties i.e. size, color, fabric etc. are available on demand. We guarantee the best value and competitive services at affordable rates.

SALIENT FEATURES

- ◆ Hospital curtain polyester fabric water repellent and stain retardant
- ◆ Hospital curtain polyester fabric fire retardant, anti microbial and anti fungal up to 50 washes
- ◆ Standard curtain size:4 x 7 feet

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(23)

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

TWO POINT IV BAG HOLDER MODEL NO.: SEP-IV01



DESCRIPTION

It is designed specifically for the overhead suspension of drip bags without the need to use valuable floor space. It is always available at a moment's notice, yet remains discreetly unobtrusive. As the carrier is track mounted, it can be moved to one side of the cubicle when not in use, and can be positioned at either side of the bed, or above it if required.

SALIENT FEATURES

- ♦ Strong aluminum track, white powder-coated finish
- ♦ Economical 2-point drip carriers with nylon runners
- ♦ Track suspended at cubicle track height
- ♦ Capacity 3 kg per hook
- ♦ Hooks made of stainless steel SS304
- ♦ Mobile track option
- ♦ Suspension system available in stainless steel S304 and MS powder coating.

TELESCOPIC IV BAG HOLDER MODEL NO.: SEP-IV02



DESCRIPTION

It is becoming increasingly important that the area around a hospital bed should be as uncluttered as possible. It is now accepted that the patient's quality of life can have a direct bearing on recovery time, and the perceived intrusiveness of medical equipment may be a contributory factor.

SALIENT FEATURES

- ♦ Heavy-Duty 5 Bag Holder
- ♦ Heavy duty self-locking wheeled trolleys
- ♦ Five point telescopic carrier trees
- ♦ Made of stainless steel SS304 and aluminum with white powder coating on it & hooks are of SS304
- ♦ Strong aluminum track with white powder - coated finish
- ♦ Track suspended at cubicle track height

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

24

© Copyright Reserved by : Western Surgical - Rajkot.



A1 2,000/-



A2 35,00/-



B1 2,700/-



B2 5,000/-



C1 4,500/-



C2 8,500/-



D1



D2



E

**ABS BED
SIDE
LOCKER**

9,500/-



F1



F2

11,000/-

WS - 24

BED SIDE LOCKER

Available in three Models.

- Locker made from machine pressed CRCA Sheet with S.S. Top having edges on three sides. Pretreated & Epoxy Power Coated/complete SS
- Locker made from machine pressed CRCA Sheet with three sides enclosed with one box and one drawer. Provided with 50 mm dia revolving castors, with S.S. Top, having edges raised on three sides. Pretreated & Epoxy Power Coated/Complete Stainless Steel.
- Trolley / made from pressed bent stainless steel Sheet with stainless steel top. (12" x 17" x 30") fully stainless steel with 50 mm dia. Revolving castors.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(25)

© Copyright Reserved by : Western Surgical - Rajkot.

WS - 31A



Cylinder Trolley Jumbo

Tubular Structure with heavy duty 200mm dia castor wheels. To support easy movement for jumbo cylinders.

2 Wheels Oxygen Trolley



Oxygen Cylinder Trolley Height Adjustable



2,000/-



WS - 31

1,200/-

Cylinder Trolley 2 wheels regular

Tubular Structure with 100mm dia rubber tyres, easy to carry 10.2 ltrs Size Oxygen cylinder. Available in complete powder coated or Stainless Steel.



Cylinder Trolley 6 wheels Regular

1,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

(26)

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

WS- 1095A

Folding Stretcher (canvas)

3,500/-



(213 L x 56 W x 15 H cms)
M.S. Pole width wise fold.

WS- 1095B

Folding Stretcher 2 Fold

5,500/-



Unfolded Size : 221L x 53W x 15 Hcms
Folded Size : 110L x 18W x 11 Hcms
Net Wt : 5 kg
Gross Wt : 6 kg
Load Baring Capacity : 159 kg

WS- 1095C

Folding Stretcher 4 Fold

7,000/-



Unfolded Size : 210L x 53W x 12 Hcms
Folded Size : 61L x 30W x 12 Hcms
Net Wt : 6 kg
Gross Wt : 7 kg
Load Baring Capacity : 159 kg

WS- 987A Scoop Stretcher (Telescopic/immobile)

A portable first aid stretcher, is commonly used for transporting bone fractured and seriously injured patients. A clutch zip design at the middle of the model, the stretcher can be divided into left and right parts. Both parts can be drawn out when you need to put the patient on the operating table or the hospital bed quickly without moving the patient. The stretcher length can be directly transferred into the ambulance. The stretcher can fold or be separated into two parts, easy to carry or transport.

Technical Parameter :

Product size (LxWxH)	: 168x43x7cm(min) : 224x45x10cm(max)
Net Weight	: 9 kg
Gross Weight	: 11 kg
The bearing pressure	: < 159 kg
Dimension of the packing box :	127x45x10cm



8,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

Name & Model

Neil Robertson Stretcher

NRS - 01

INR 10,875.00 + 2% Against C-Form or 12.5%.

Delivery within 6 weeks after the order is received.

Material & Specification

High tenacity PVC fabric and wood



12,500/-

Ambulance Stretcher

MS Top & S.S. Top with Bearing wheel as per your Ambulance



3,500/-



8,500/-

WS - 34A

FOLDING TYPE STRETCHER TROLLEY

Overall Size : 72" x 22" x 32"

Tubular steel structure with both sides folding on accessible rack mechanism. Optionally available with 4" or 7" dia. Wheels. Pretreated & Epoxy Power Coated or completely Stainless Steel.




4,500/-

Folding Stretcher New Design



9,500/-



Stretcher Trolley Folding With Oxygen Stand

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

(28)

© Copyright Reserved by : Western Surgical - Rajkot.

WS - 34**STRETCHER ON TROLLEY****Overall Size : 84" x 22" x 32"**

Tubular steel structure with removable stretcher of M.S. Sheet Top. Optionally available with 4" or 7" dia. Wheels. Oxygen Cylinder Cage, Tray for Instruments, I.V. attachment with saline Rod. Pretreated & Epoxy Power Coated or completely Stainless Steel.

M.S. TOP :- 6,000/-
S.S. TOP :- 7,500/-

WS - 33A**SPECIAL STRETCHER ON TROLLEY****Overall Size : 84" x 22" x 32"**

Tubular steel structure with removable stretcher of M.S. Sheet Top. available with 4" Whels. Oxygen Cylinder Coge, Tray for Instruments, I.V. attachment with saline Rod. Pretreated & Epoxy Power Coated



15,000/-

WS-34B**Stretcher on Trolley****WS-814****Patient Trolley (Hi-Low)**

- ◆ Approx. size : 195L x 65W x 68-90H cms. ◆ Frame Work made of S. S. Tubes
- ◆ S. S. Side railing, I.V. rod. ◆ Adjustable Height by crank handle ◆ 15 cms dia castors.
- ◆ Full Stainless Steel. ◆ Also Available in complete M.S. frame.



25,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.




WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(29)

© Copyright Reserved by : Western Surgical - Rajkot.



21,500/-

WS - 32

TRANSFER TROLLEY SYSTEM

"The exchange area technique has come to stay. This new concept prevents any possibility of cross infection. On the way to the operating suite the patient is brought to the exchange area on a ward trolley. The under carriage of the ward trolley and that of the operating room along together in the exchange area. With the trolley top the patient is transferred to the operating room undercarriage on easy rolling castors. The undercarriage is then separated and patient is taken to the operating room. The reverse process taken place after the operation when the patient is to be shifted to his ward.

Tubular steel structure with removable stretcher of M.S. Sheet Top. Available with 7" dia. Wheels. Oxygen Cylinder Cage, Tray for Instruments, I.V. Attachment with saline Rod. Pretreated & Epoxy Power Coated or completely Stainless Steel.

WS - 1008

Ambulance Stretcher

Light weight, strong structure, lovely design, good arts and crafts, easy antiseptic, small size, convenient to use. Utilizes high precision mechanical transmission. Can be used automatically folded through the left and right hand control handle. Made of high quality aluminum alloy tubes. Dimensions :



33,000/-

Stretcher Trolley CT & MRI



45,000/-

A - High position : 1,900 x 560 x 900mm
 B - low position : 1,900 x 560 x 280mm
 Load weight : upto 120Kg
 Product weight : 40Kg
 Diameter of the castors : 125mm
 Applicable for various kind of hospitals, emergency centers, ambulances and battle fields

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54

30

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

Patient Transfer Trolleys

9,500/-

USI- 991

Patient transfer system

Patient mover from bed to bed

Applicable in large, medium and small hospitals, ICU, Operation Theatre, Emergency Room and Patient Ward, etc.

Usage:

According to current requirement of modern nursing science, nursing involves a lot of physical work, which causes backache and lumber pain on parts of the nurses. While being used to move the patient, especially the patient who needs special care, between stretchers, patient beds, trolleys inspection beds, X-ray fluoroscopy tables and operation tables, EASYMOVE is easy to use.



Technical Specifications:

Description	: Patient mover from bed to bed
Size	: 170 x 50 cm
Material	: Nylon, Plywood, Water-resistant.
Features	: Easy to be clean and sterilized. Maximum load for 15cm gap will be 200kg.



NURSES WORKS SIMPLER
PATIENTS FEELS SAFE

WS - 985 Wheelchair Stretcher with varied positions

35,000/-

Technical parameter :

(Product size) (LxWxH): 190x59x85 cm

(Net weight): 45kg

(The bearing pressure): <159kg.

(Dimension of the packing box):

200x60x27cm

(Gross weight): 55kg



Automatic wheelchair type stretcher is mainly used to transport the sick, wooded outdoor. it's a necessary accessory on ambulance and could be folded as wheelchair in narrow stairs or elevators. This model is made of high quality aluminum alloy, with characteristics such as light weight, good appearances, reasonable configuration, convenience operation, antiseptis. Besides, it has soft foam matters block sidely. it's an easy operation stretcher.

WS - 34Z

TRAUMA CARE

RECOVERY TROLLEY

Overall Size : 84" x 24" x 24-32"

- Height : Adjustable by foot operated Hydraulic pump
- X-Ray Permeable removable stretcher with built-in back rest
- Trendelenburg / Reverse trendelenburg position on gas spring
- Swing away Stainless steel side rails
- Oxygen Cylinder Holder
- Sliding X-Ray Cassette tray



65,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

31

© Copyright Reserved by : Western Surgical - Rajkot.

WS- 988 Spine Board

9,500/-



Spine Board is mainly used for immobilizing and safely transferring serious patients during rescuing. It's a first-aid transport equipment for the medical department. It also could go into X-ray directly and float on water. This model is made of high-strength engineering plastics, it can be enduring used, and has the advantage of anti-aging, weight loading, and easy stocking specialties.

The Features of our Spine Board are :

- CE certified Spine Board- in line with European Standard
- Withstands up to 300 kgs.
- Available with 3 belts & groove to hold the board
- Rigidity-Long life and light
- Light to carry
- Floatation
- Resistance to bump and corrosion
- Does not absorb Liquids
- Immune to infiltration
- Easy to clean (Water and soap are enough)
- Translucent
- Homogeneity with X-Ray
- Seamless moulding to eliminate ingress of fluid
- Avoids splintering and cracking
- Rigid Support During Movement of a Patient with Suspected Spinal or Limb Injuries., Ambulance Services and Lifeguards
- Better in Strength, User Friendly & Comply with Exacting Global Medical Standards

Specification of Spine Board

Sr. No.	Particulars	Standard	Test Method
1	Description	HMHDPE plastic molded stretcher for immobilization of patient.	Visual
2	Raw Material	Virgin HMHDPE granules form Lotrene Having grade TR-571 for Spine Board and Nylon for Belt.	Test Certificates
3	Colour	Orange	Visual
4	Total Article Weight (Kgs)	7.200 +/-2%	Electronic Weighing Balance
	Dimensional		
1	Article Length (mm)	1823 +/-5mm	Measuring Tape / Scale
2	Article Width (mm) Max.	448 +/-5mm	Measuring Tape / Scale
3	Article Thickness (mm)	61 +/-5mm	Measuring Tape / Scale
4	Load Capacity (Kgs.)	300.00	TM/QC/02

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



© Copyright Reserved by : Western Surgical - Rajkot.

WS - 42**WASTE CARRYING TROLLEY**

Trolley to carry waste basket made in complete plastic. Trolley completely in Stainless Steel with 100mm dia heavy duty revolving castors. M.S. :- 3,0



5,500/-



**MOPING
TROLLEY
3 Dustbin**

WS-808**Bio Medical Waste Bins**

Available in 100 Ltr. Capacity with stand without foot operated.




Size Available : 12 Ltr. / 30 Ltr. / 80 Ltr.

6,000/- 7,500/- 10,500/-

Mopping Trolley 3 Bucket Systems**Bio Medical Waste Management & Handling**

13,500/-

BMW Box Type Trolley**MS / Steel Both****Available Size :****L 48" x H 30" x W 30"**

35,000/-



38,000/-

**CSSD Air
Tight Trolley**

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

33

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



6,500/-

Reg.:..... Del.:.....

WS - 14

EXAMINATION TABLE HEAD UP

Overall Size : 72" x 24" x 30"

72" x 22" x 30"

Table top in two section. Head Flap adjustable on a ratchet. Pretreated and powdercoated Optionally available Stainless Steel Top or completely Stainless Steel.

WS-EXAMINATION TABLE WITH DRAWER

Reg. :.....Del.:.....



8,800/-

WS-EXAMINATION TABLE WITHOUT DRAWER

Reg. :.....Del.:.....



6,000/-

WS-EXAMINATION COUCH WITH HEAD UP

19,500/-



A

Examination Table (Electric)

Model No. : WS-1060

Overall Size :

1800L x 600W mm

Dimension (Head)

560L x 600W mm



WSGEC2D3LM : Gynec Examination Couch 2-D-3LM



24,500/-

WS-141 EXAMINATION CUCH GYNAEC TYPE



35,000/-

PHYSICIAN EXAMINATION CUCH



20,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

34

© Copyright Reserved by : Western Surgical - Rajkot.

S.S. TOP :- 15,000/-



WS - 15

OPERATION TABLE
SCREW TYPE

Overall Size : 72" x 22" x 32"

Table top in three sections provides trendelenberg, reverse trendelenberg, which can be obtain with a simple screw mechanism. Leg flap can be detached, provided with lithotomy rods, handrest and single step stool. Pretreated & Epoxy Powder Coated. Completely Stainless steel

WS - DELIVERY TABLE STEEL TOP

Size : 72"x30"x30" ht.



11,500/-

WS - DELIVERY TABLE MS TOP

Size : 72"x30"x30" ht.



DELIVERY TABLE FULLY STEEL



WSGEC12-DM :

Gynec Examination Couch-12D



30,000/-

DELIVERY TABLE



WS 142 :

EXAMINATION TABLE GYNAEC.

Size : 1830 x 560 x 760 mm ht.

72" x 22" x 30" ht.

Available in

- a) MS Powder Coated
- b) S. S. 304 grade.



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

35

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

WS - 203 TELESCOPIC DELIVERY TABLE

Size :72"x 30" x 30" ht.



Tubular frame work powder coated or fully S.S. with 18 Swg. stainless steel Top. Provided with saline rod, thigh rest, stainless steel bowl. Telescopic foot end on smooth heavy duty sliding channels, Backrest adjustable on ratchet and head low with gas spring mechanism. Provided with powder coated or S.S. Detachable side rails and head low.

Attendance Bed with Locker Table



6000/-

WSAB2L : Attendance Bed & 2 Locker




7600/-

WSAB-4LWM : Attendance Bed - 4L Without Mattress



8500/-

WSAB4LM : Attendance Bed & Mattress



10,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

36

© Copyright Reserved by : Western Surgical - Rajkot.

Delivery Table

Gynceologic Couch

Electric



1. Overall Size : 1900mm L x 620mm B x 700-950mm H .
2. Frame work made up of CRCA tabular pipe frame.
3. Three section cushioned mattress top with detachable leg section.
4. Trendlenburg / Reverse Trendlenburg Positions , Height Adjustment, Backsection adjustment through Actuator System.
5. Provided with a pair of Upholstery Lithotomy Positions
6. Detachable S.S Tray.

Remot :- 85,000/-
Hydraulic :- 65,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

(37)

© Copyright Reserved by : Western Surgical - Rajkot.

Delivery Table Hydraulic



90,000/-

1. Overall Size : 1900mm L x 900mm W x 600-800mm H
2. Tabular Frame Work mounted on 125mm dia castor two with brakes.
3. Three Section top made of M.S Pipe with fixed cushioned mattress.
4. Trendlenburg / Reverse Trendlenburg Positions adjustment through Gear Mechanism.

5. Height Adjustment through Imported Hydraulic Pump.
6. Back section Adjustment through ratchet system.
7. Leg end section can slide under the main section.
8. Provided with a pair of Upholstery Lithonomy Positions
9. S.S Telescopic I.V Rod.
10. Detachable S.S Tray.



Electric

1. Overall Size : 1880mm L x 900 mm W x 600-800mm H
2. Tabular Frame Work mounted on 125mm dia castor two with brakes.
3. Three Section top made of M.S Pipe with fixed cushioned mattress.
4. Trendlenburg / Reverse Trendlenburg Positions , Height Adjustment, Backsection Adjustment through electric Actuator System.
5. Leg end section can slide under the main section.
6. Provided with a pair of Upholstery Lithonomy Positions
7. S.S Telescopic I.V Rod.
8. Detachable S.S Tray.

1,50,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

38

© Copyright Reserved by : Western Surgical - Rajkot.



WS - 45

INSTRUMENT TROLLEY S.S. TOP

Overall Size :	27" x 18" x 30"	30" x 20" x 30"
	4,500/-	6,000/-

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS




WS - 47

DRESSING TROLLEY STEEL

Overall Size : 30" x 20" x 30"

Tubular structure with two stainless trays, having side rails with brackets for stainless bowl and bucket. Mounted on 100 mm dia. Revolving castors. Pretreated & Epoxy Power Coated or fully S.S.

9,000/-



INSTRUMENTS TROLLEY FULLY S.S. STEEL

Overall Size :	27" x 18" x 30"	30" x 20" x 30"
	7,500/-	9,500/-

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



**WS - 245
RECT. STAINLESS
STEEL BED
SIDE TROLLEY**

Overall Size
430 x 300 x 760mm ht.

6,500/-



**WS - 246
SQ. STAINLESS
STEEL BED
SIDE TROLLEY**

Overall Size
380 x 380 x 760mm ht.

6,800/-



21,000/-

Cathlab Trolley SS

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



INSTRUMENT TROLLEY S.S.

Overall Size : **27" x 18" x 30"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



INSTRUMENT TROLLEY S.S.

Overall Size : **30" x 20" x 30"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



INSTRUMENT TROLLEY S.S.

Overall Size : **48" x 24"**

Tubular Structure Mounted on 50 mm dia. Revolving castor with stainless steel top and shelf. Pretreated & Epoxy Powder Coated.

Also Available in MS / SS



3 SELF STEEL TROLLEY

Overall Size : **20" x 30"**

A-Steel Top
B-Fully Steel



3 SELF STEEL TROLLEY

Overall Size : **18" x 15"**

A-Steel Top
B-Fully Steel

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Any Requirements Sizes Are Also Available We Manufacturing



7,200/-

ECG TROLLEY STEEL

18"x18"x33"h
fully steel 3 drawer Modolar

FOOD TROLLEY STEEL & MS 24" x 48" x 2 Self



25,000/-



"Z" Type Trolley
3 Self MS &
S. S. Top &
Fully Steel

M.S. 4,500/-
S.S. 7,200/-

ECG / DIFEBRILATOR TROLLEY



3,500/-

E.C.G. TROLLEY WITH TWO SELF

DERMA TROLLEY



12,500/-

Defibrillator Trolley



7,200/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



7,500/-

9,500/-

WS - 71

BABY RESUSCITATION TROLLEY

Overall Size : 27" x 18" x 30"
30" x 20" x 30"

Stainless Steel tray in head low position.
Three drawer modular system.
Adjustable I.V. Rod, Oxygen cylinder cage & lamp arrangement.
Bottom tray for mounting suction apparatus.
Trolley mounted on 100 mm dia wheels.
Pre-treated & epoxy powder coated.
Also available in stainless steel.



6,500/-

**ECG Trolley
Steel Top 3 Self with
Drawer
18"x18"x33"h**



**ECG Trolley
Fully Steel 3 Self
18"x20"x36"h**



7,200/-

**ECG Trolley
Fully Steel 3 Self with
18"x20"x36"h**



5,800/-

**ECG Trolley
Steel Top 3 Self
18"x18"x33"h**



**ECG Trolley
Fully MS 3 Self
18"x18"x33"h**



**ECG Trolley
Fully SS 3 Self
18"x18"x33"h**



**ECG Trolley
Fully SS 3 Self
18"x18"x33"h**



**Automatic X-Ray
Processor Trolley
18,000/-**

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



9,500/-

EMERGENCY TROLLEY WITH DRAWER



9,500/-

EMERGENCY TROLLEY WITH DRAWER

Laundry Trolley Required Specification



Size : 660 x 1600 x 1676 mm Stainless Steel, 2 Shelves
Moves Easily on Four 6" Heavy Duty Casters

RETINA MACHINE TROLLEY



15,500/-

Size : 20" x 27" x 33" x 4
Self with 7 Electric Point Panel



WS - 131

FRESH & SOILED LENIN TROLLEY
Size : 460 x 1090 x 890mm ht.



WS - 132

SOILED LENIN TROLLEY EXCELL
Size : 790 x 560 x 1030mm ht.



WS 132 : SOILED LINEN TROLLEY EXCEL.
Size : 790 x 560 x 890mm ht/ 31" x 22" x 35"
100mm / 4" non rusting, non marking castors.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

43

© Copyright Reserved by : Western Surgical - Rajkot.

Mayos Hight Adgestable Trolley Complete Steel

Size : 20" x 30" x 30" / 45" H/A
Also Available Costumer Size



Lasar Machine Trolley For Ophthalmic



27,500/-

V.H. Trolley S.S. Top & Comp. S.S.

Size : 18" x 18" x 20"



Mayos Trolley

Complete Steel
Size : 18" x 21" x 30" / 45" H/A
Also Available Costumer Size



6,000/-

WS - 40

Mayo's table

Overall Size : 22"x 16" 32" - 50"

Tubular Steel Structure with Stainless Steel deep Tray. Height adjustable with the help of telescopic S.S. Tube. Pretreated & Epoxy Power Coated.

Also available in fully S.S.

3,800/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

WS-38 OVER BED TABLE PUSH TYPE

44" x 17.5" x 38" ht

Tubular Steel Structure with laminated Tray on 50mm dia, revolving castors pre-treated & powder coated



4,100/-

WS - 36

ADJUSTABLE BED SIDE TABLE

7,200/-

Overall Size : 24" x 18" x 32" - 50"

Height adjustable, tubular steel structure with laminated top, Adjustable on S.S. Telescopic tube. Mounted on 50mm dia. Revolving castors. Pretreated & Epoxy Powder Coated. Also Available in Fully Stainless Steel. Optional available fully s.s.



ADJUSTABLE BED SIDE TABLE GEAR TYPE

Overall Size : 30" x 16" x 32" - 50"

Height adjustable with help of gear mechanism provided on the side of the table. Tubular steel structure with laminated top with adjustable tilt. Mounted on 50mm dia. Revolving Castors. Pre treated & Epoxy Powder Coated. Also Available in Fully Stainless Steel.

7,200/-

WS - 36 A



Available in Attractive Colours



WS 601 OVERBED TABLE FOLDING WITH HYDRAULIC FUNCTION



6,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Scrub Station & Instrument Washing



Dimension's As Per Your Size



Dimension's As Per Your Size



4 Feet

Dimension's As Per Your Size



65,000/-

WS - X13

SCRUB SINK

Size : 20" d X 39" h X vertical depth 18"

Scrub sink made in complete 18 swg 304 grend SS Sheet. Tailor made to fit any lenght, minimum 30".

Drainage options either on left or right Available with or without bottom cover Can be provided with foot pedal taps automatic sensor taps. Soap Dispensers



Dimension's As Per Your Size

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



3,600/-
Fully Steel



1,800/-
S. S. Top



1,500/-
S. S. Top



1,500/-

WS - 39

REVOLVING STOOL

Overall Size : 11.5" dia x 21" - 29"

Tubular Steel Structure with Superimposed Stainless Steel Top. Height adjustable from 425 to 675 mm with the help of smooth screw mechanism. Pretreated & Epoxy Power Coated. Optionally available fully Stainless Steel.



4,100/-
Hydrolic with Wheel



3,000/-
Hydrolic Fix



S. S. Top
1,500/-



M.S. Top
1,200/-



Fully Steel
2,700/-



WS 393 : REVOLVING STOOL WITH CUSHION.
Size : 370 dia 545-745mm ht. / 14" dia 17"-29"
Available in MS Powder coated / S. S 304 grade.

WS - 68

BED SIDE STOOL (Overall Size : 12" x 12" x 8")

Mode from CRCA steel structure with M.S. Pretreated Epoxy Power Coated. Optionally available with laminated top, Stainless Steel top or fully Stainless Steel.



3,200/-

Food Trolley on Bed
Size : 18" x 24" x 15" H.
Complete Steel & Powder Coted



WS 393 : REVOLVING STOOL WITH CUSHION.
Size : 370 dia 545-745mm ht. / 14" dia 17"-29"
Available in MS Powder coated / S. S 304 grade.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

WS - 13

SOILED
LENIN BIN

Overall Size : 20" dia X 36" ht

Tubular Structure mounted on 100mm dia revolving castors. Provided with detachable canvas back. Pre-treated epoxy powder coated / Complete Stainless Steel.

2,500/-



WS - 15 G

OPERATION TABLE FOR
ANIMAL SURGERY

HEIGHT ADJUSTABLE

Size : 48" x 26" x 30" ht.

Epoxy Powder Coated. Completely Stainless steel

11,000/-



KICK BUCKET SMALL

2,100/-



S - 54

WS - 57

BOWL STAND

Overall Size : 33" X 14" dia with

Tubular Structure mounted on 50mm dia. revolving castor with two stainless steel bowl. Pretreated & Epoxy Powder Coated / completely SS.



2,700/-

WASH BASIN SINGLE

Overall Size : 35" X 14" dia with

Tubular Structure mounted on 50mm dia. revolving castor with a stainless steel bowl. Pretreated & Epoxy Powder Coated / completely SS.



2,600/-



3,600/-

**Wash Basin
MS & SS**

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *





Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

48

© Copyright Reserved by : Western Surgical - Rajkot.



**WS W1
WAITING CHAIR
EXCEL CUSHION.**

- W 1-1 Single Seat with PU Cushion
- W 1-2 Double Seater with PU Cushion
- W 1-3 Three Seater with PU Cushion
- W 1-4 Four Seater with PU Cushion

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54



15,000/-

WS - 66

SURGICAL INSTRUMENT CABINET

Size : 78" X 19.5" d X 36" W

Compact Cabinet for instruments provided with 4 shelves to form 5 compartments. Front doors in glass 5mm thick.



18,000/-

WS - 67

INSTRUMENT CABINET

Size : 63" X 17.5" d X 25" W

Frame : CRCA rectangular tube, CRCA sheet plants
Five Shelves, door & two sides of 5mm thick fine glass. Pre-treated & epoxy powder coated.

Surgical Drawer

Surgical Cupboard

Multi Locker - 4

Multi Locker - 12



Multi Locker - 12

WS 673 : CABINET 6 COM.

WS 674 : CABINET 15 COM.



Size : 1970 x 500 x 920mm ht.
78" x 19" x 36"
0.8mm thick CRCA sheet
Completely Powder Coated.



Size : 1970 x 500 x 920mm ht.
78" x 19" x 36"
0.8mm thick CRCA sheet
Completely Powder Coated.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Bed Side Screen - 3 Fold



3,500/-

WS - 67 BED SIDE SCREEN

Size : H 66" X W 96"

Middle Span of 48" tubular structure in 3 sections mounted on 50mm revolving castors. Available in complete Stainless Steel / Power Coated.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

51

© Copyright Reserved by : Western Surgical - Rajkot.



WS-

4,500/-

WS-

5,000/-

4,500/-

6,500/-

Surgeon Stool

Tubular Steel Structure with Superimposed Stainless Steel Top. Height adjustable with the help of Heavy Duty gas spring. Back made up of uniformly perforated CRCA. Sheet. Pretreated & Epoxy Power Coated.



WS- 1055 Surgeon Chair

Overall Approx. Size : Cushioned seat and back
Pneumatic height adjustable
Height adjustable and swinging back rest
5 Legs chrome plated base with foot ring
Base fitted with 50 mm dia castors

WS - 59

Overall Size :
20"x18"x18"x20"x12"x8"
Tubular Structure with M.S.
Top with Aluminum
Treads. Pretreated
and Epoxy Powder
Coated.



Single Step
Stool M.S.



Double Step
Stool M.S.



Single Step
Stool Fully Steel



Double Step
Stool Fully
Steel

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Injection Trolley 3 Self + 3L



12,500/-



27,000/-

Crash Cart Steel Trolley



Crash Cart Trolley Mini with 12 Drawer



18,500/-



20,500/-



24,000/-



Emergency I.C.U., Defibrillator, Monitor & Suction Trolley
Size : 27" x 27" H 48" 14,000/-

WS - 69 EMERGENCY CRASH CART

Tubular trolley with two boxes Alkon make with three drawers each, removable multicolored plastic bins with saline attachment and cylinder trolley. Provided with 100 mm dia. Revolving Castors, Pretreated and Epoxy Powder coated.

Optionally available fully stainless steel, cardiac message board & light source.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54



© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

Emergency I.C.U., Defibrillator Trolley

21,000/-



No. 1



No. 3

Size : 27" x 27" Two Drawer with Partition

12,500/-



35,000/-

WS - 473 TRAUMA CARE MEDICINE TROLLEY

WS - N7 EMERGENCY DRUG TROLLEY



A innovative product with multicolored compartments in each row provided with additional 3 drawers on heavy duty telescopic channels. Handling made easy with Handle provided at the rear. A big compartment assessable from behind os Useful to store bulk Material.

MODELS : A) 30 Compartments / B) 50 Compartments
(A) 23,000/- (B) 35,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

54

© Copyright Reserved by : Western Surgical - Rajkot.

Endo / Laparoscopy Trolley



Special Features

1. Size - 4.75 Ft. H x L + W - 2.25 Feet
2. Attractive Color Combination
3. Self Size - 18" x 24" tray type
4. Self Height Adjustable No. 2,3,4,5
5. One Drawer with Lock
6. Stand for Cauty foot switch
7. 4" Castor with Lock systems
8. Two S.S. Handle for Movement
9. Co2 Cylinder stand Back Side with 2 clamp
10. 6 + 1 Plug 15 x 5 Amp Switch & one Extra for Cauty with 10 Ft. Cable
11. Space for UPS / Servo Stabilizer behind



DUAL CABLE CONNECTION

18,000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

55

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

Economy Laparoscopy / Monitor Trolley



Special Features

1. Size - 4.75 Ft. H x L + W - 2.25 Feet
2. Attractive Color Combination
3. Self Size - 18" x 20" tray type
4. Self Height Adjustable No. 2,3,4
5. One Drawer with Lock
6. 4" Castor with Lock systems
7. Two M. S. Handle for Movement
8. Co2 Cylinder stand Back Side with 2 clamp
9. 6 Plug 15 x 5 Amp Switch with 10 Ft. Cable
10. Space for UPS / Servo Stabilizer behind



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

56

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



5,500/-

Perforated E.T.O. Tray For C.S.S.D.
Size : 18" x 15" / 18 G / 20 G / 22 G



6,500/-

Net Type E.T.O. Tray For C.S.S.D.
Size : 18" x 15"



5,500/-
 15,000/-

Cidex Tray for Lapro/Urology

Size : 1) 27"x6"x5" (2 Compartment) - Lapro 2) 44"x9"x6" (2 Compartment) - Cardiology
3) 32"x8.5"x4.5" (2 Compartment) - Uro/Ortho 4) Steel / P.V.C. / Acrylic Available Tray



9,500/-
 35,000/-

Cidex Tray 2 in 1 Lapro/Urology

Size : 1) 27"x6"x5" (2 Compartment) - Lapro 2) 32"x8.5"x4.5" (2 Compartment) - Uro/Ortho
3) 44"x9"x6" (2 Compartment) - Cardiology 4) All Steel Tray


* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54

BURN CAGE



3,500/-

Catheter Hook For Cardiac OT-Lab



3,500/-

Variable Height Tilt Table



Height adjustable tilt table with two motors. One motor ensures effortless and smooth tilting; the other adjusts height and facilitates patient transfer. The Variable Height Tilt Table is the ideal standing aid for neuromotor revalidation and adaptation exercises from lying to standing.

Another area where it can be applied is general orthostatic examination, more specifically examination and treatment of vasovagal syncope. Standard with three straps.

Standard model comes with a dual adjustable footplate. The tilt table also has a hoist access as standard and a quick release Trendelenburg positioning. It can also be fitted with an optional Torso Support for improved patient support and comfort, as well as an activity tray - for extra working space for hand activities.

Dimensions



Height



Table Positions



Specifications

- Minimum Height : 51cm/20"
- Maximum Height : 86cm/33.9"
- Cushion Width : 65cm/25.6"
- Leg Length (+ Footplate) : 200cm/78.7"
- Lift Mechanism : Electric
- Breathing Hole : 0 / Wheels : 4 Braked
- Tilting Mechanism : Electric
- User Weight Rating : 205kg/451lb



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

58

© Copyright Reserved by : Western Surgical - Rajkot.

WS - N4

DRUG TROLLEY TWIN

Size : 18" X 20" X 48"

Most compact drug trolley having 3 different sizes of heavy duty bins placed louvered panel making it adjustable. 23 bins on each side. Additional tray provided at bottom. Can be provided in single side also.



24,500/-

WS - X15

X-RAY APPRON STAND



11,000/-

Syringe Pump Stand



Patient Shifting Rollar



5,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

59

© Copyright Reserved by : Western Surgical - Rajkot.

WS- 939

Emergency Trolley (ABS)

ON DEMAND

ABS body, metal ball-bearing rail. A-Aluminium column B-Plastic steel columns. IV Pole, needle disposal holder (optional); Dust basket, utility container, file bag, defibrillator shelf, power outlet & hooks, Oxygen tank holder (Optional), C.P.R. Board (Optional), Sliding side shelf, stainless steel guard rail, concave top covered transparent soft plastic glass, Centralized lock, 2 small & 3 big drawers repositioned easily with partitions can be well organized. ABS basket & drawer frame can be repositioned easily. Luxurious wise casters, two with brakes.

Approximate Size : 850 x 520 x 920 mm



WS- 940

Medicine Trolley (Cassette Type)

- * Designed to supply medicine when nurse make around of visit to patient
- * Two sides unit cassette and one side plate cassette
- * Easy to discharge
- * Installed multi-bin container at rear side to make easy and multi possible strong
- * Multi structure can be possible used by kinds of accessories
- * Ergonomics designed push-handle and efficient structure
- * Total lock key system

Size : 780(L) x 490 (D) x 920(H)



ON DEMAND

WS- 40D

Case History Trolley

Approx Size : 680 x 400 x 850 mm

Made of ABS Material

Double row, 40 layers with a drawer

3" Castor

ON DEMAND



WS- 20D

Case History Trolley

Approx Size : 336 x 400 x 850 mm

Made of ABS Material

Single row, 20 Layers with a drawer

3" Castor



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54



© Copyright Reserved by : Western Surgical - Rajkot.

Blood Collection Pathology Chair



WS - S7

BLOOD TRANSFUSION CHAIR



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

61

© Copyright Reserved by : Western Surgical - Rajkot.

Overall Size: - 1830mm L X 610 mm W X Hgt. 380 mm(Approx.)

Construction: Bed /Chair frame is made up of double 25 mm dia x 18G CRCA tubes Welded to each other so cross section looks like a 50 x 25 mm tube

The top to made with built in cushion 50mm thick P.U foam and Rexene cover fixed on a ply 12 mm thick and secured to bed with help of wooded screws. The mattress should have built in pillow.

Legs of the chair should be made from 1.25" dia x 18gauge ERW tube with horizontal tube 25.4 mm x 18 swg tube Legs Provided with Nylon shoes having metal inserts The bed should easily fold and slide to make a chair and vice versa

The base to be provided with 4 fixed direction castors underneath the first section, whereas the back section of the seat should have 70 mm dia castors with brakes.

The folding section should have high gas spring to secure the chair in open position.

Addition support at the extreme front end to provide extra support and safety to any sitting at the extreme end of the bed. The hinges should be rotatable within the frame and not external.

Finish: All components should be pretreated in separate eight tank process for better finish, good adhesion and corrosion protection. Process includes Hot Degreasing, Derusting, Activation, Phosphating& No's of Water rinses as per IS 3618 - 1966 class 'C' type and then pretreated materials is coated with epoxy powder with film thickness of 60 microns (approx.) and then oven baked at 180 degree centigrade.

ATTENDENT BED / CHAIR



9,000/-



WS 604 : SOFA CHAIR TWIN SEATER

Overall Size :

1980 x 1035 -810mm ht.

78" x 40"x 32" ht.

Mattress : 1980 x 880mm / 78" x 34"



WS 605 : SOFA CHAIR SINGLE SEATER

Overall Size :

800mm x 710mm x 430mm ht.

32" x 28" x 17"

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

62

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

QUADRICEPS EXERCISE TABLE



7,500/-



15,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159


Mo. +91 94280 03750/52/54

63

© Copyright Reserved by : Western Surgical - Rajkot.




D1 - 360 x 275 x 17




Central locking system is provided OPTIONALLY.




D2 - 360 x 275 x 40



Special Stoppers are provided to stop the Drawers when open fully



D3 - 360 x 275 x 63





To change the name plate, push the indicator in upward direction.



D4 - 360 x 275 x 88




MIXED DRAWER SYSTEM
Modular drawer system can be used with different sizes of drawer in the same cabinet.



ART No.	A	B	C
SB-1	125	100	60
SB-2	170	100	85
SB-3	220	130	125
SB-4	290	178	143
SB-5	290	363	145
SB-6	440	180	145
SB-7	435	400	215

Rear lip for easy access.
Stops bin from coming out.



Art No. ACO 25



ACO 15/25

Special Holes Are Provided On The Backside Of The Unit For Wall Mounting.



Special Holes Are Provided On The Backside Of The Unit For Wall Mounting.

Available in Two Colours :

DRAWERS Transparent	CABINET Grey colour
-------------------------------	-------------------------------



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Shoe Rack Steel



All Size Different Available

Shoe Rack Steel 48" x 12" x 36"



All Size Different Available

Medical Rack For Surgical Medicine File



Medical Rack For Surgical Medicine File



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

65

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

MEDICAL STUDENT CHAIR



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

66

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

**Medical
Equipments
&
Instruments
Safety
Locker**



**fire and Burglar
Resisting safe,
Deposit safe Cash
Box, Coffe, Strong
Room Doors, Safe
Deposit Lockers,
Safe Cum Safe
Deposit Locker,
Fire Resisting
Record Cabinet,
Fire Resisting Filing
Cabinet, Fire
Doors etc.**



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED**

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

Western Surgical

DOCTOR'S CHAIR



Doctor Office Table

7,500/-

5,500/-



**Medical
Computer
Table**



4,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

68

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

HOSPITAL VISITOR CHAIR



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

69

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

Customized Plastic Moulded Chair Parts



- ❖ Impressive aesthetics ❖ Robust designs
- ❖ Global quality standard ❖ Durable ❖ Smooth finish

Customized Molding Products



Plastic Chair Parts
(C INFC)

Plastic baby car
seat shell (C BSS)

BUS Seat
(AC-10)

Plastic stadium seat
(C Stoll)

Evaluating capabilities and strengths coupled with market needs led to the origin of this niche business unit, Infrastructural Furniture: with finish and quality at par with global standards. Already a product that is well accepted by leading seating systems manufactures; the expansive range includes blow moulded stadium seats, seat shells for bus, seat and back rest for educational furniture & cafes and bus seats etc.

Standard



PLASTIC CHAIR SHELLS (C INLP)

FEATURES :

- Smooth Finish at par with global standards
- Easy fitting to the frame
- Available in different colors on demand
- Aesthetically appealing


* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

CURTAIN TRACKS




DESCRIPTION

Based on the talent and competency of our skilled workers, we bring forward a distinctive collection of Hospital Curtain Tracks. In order to attain maximum client satisfaction, we are providing these tracks in various size and dimensions. Used for hanging curtains, these tracks are appreciated for their smooth and tough built. Our offered range is widely appreciated in the markets for their glossy finish, stylish pattern and vibrant color options. In addition to this, the offered range is fabricated and designed using cutting-edge techniques and superb quality fabric that add a classy look.

SALIENT FEATURES

- Durable and longer working life
- Various dimensions
- Simple to connect and operate
- Eye-catching look



CURTAIN ACCESSORIES

1. Top Bottom Bush: It is used to hang curtain tracks from the roof and made of aluminum with track color powder coating on it.
2. Curtain Runner: Made of plastic with SS hook and plastic hook and full plastic runners.
3. Curtain Removal Point: It helps to remove the curtain runner. They are fixed in track and made of mild steel with zinc plating on it to prevent rusting.
4. Jointer Plates: It helps in joining two curtain tracks together. It is made of aluminum with powder coating on it.
5. Tracks End Cap: Curtain track and caps fully made of nylon material.

VERTICAL WINDOW BLINDS

DESCRIPTION

We are offering our clients an exclusive assortment of Vertical Blinds which is extensively demanded in homes, offices and other structures. Apart from enhancing the beauty of the place, these interior Vertical Blinds prevent moisture and dust from entering into rooms. Being a customer centric organization, we offer these sturdy and low maintenance Interior Vertical Blinds at competitive prices.



BLACKOUT ROLLER BLINDS

We offer an extensive array of Blackout Roller Blinds that is available in different options and accessories. These are designed uniquely by our efficient professionals to assure flexibility in window treatment. Further, these are acknowledged for their distinctive features like flame retardant finish and user-friendly interface.

UPVC DOORS AND WINDOWS

We offer a wide variety of UPVC products including casement doors and windows, fixed windows, tilt and turn windows, French doors and windows, top hung windows, etc. The uPVC windows are sound proof, eco-friendly, insulating, termite free and low maintenance.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

POWER COLUMN

Power Column samples the organization of medical service equipment in one aesthetically appealing vertical headwall. It is a modular designed Product to provide a clean, sanitary and ordered environment. Furthermore it can enhance the quality in care and medical environment. Due to its systematized design, a network provided for medical care service is applicable to the use of secondary instruments and assemblies. The main structure consists of fluorescent panel, aluminum component and side track. Also includes power lines and medical gas piping. Side tracks are installed at four at angles, and the length is equal to the column. All the assemblies required by medical care personnel can be easily and securely placed on the side track of this column. Applications in mid to acute level environment



Features

- ◆ Designed for ease of maintenance and infection control.
- ◆ Built as a four-sided unit that offers 360° access to services and equipment management for critical care areas.
- ◆ Can be customized and upgraded to meet your facility's needs.
- ◆ Offered in a variety of widths and heights, tailored to meet your space needs.
- ◆ Wide selection of laminate colors are available.
- ◆ Multiple rail system allows for superior accessory management.
- ◆ Removable laminate panels allow for simplified maintenance.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

72

© Copyright Reserved by : Western Surgical - Rajkot.

WSP-108-A
Surgical / Anaesthesia
Pendant



Surgical Pendant Is New Addition In Field Medical Equipment In Modern Operating Theater : This Medical Pendant Has A Compact Structure By Logical Design. The Arm Console Head Can Be Rotated 340° In Low Location The Medical Staff Can Operate All The Medical Equipment In The Console Head And On The Medical Equipment Shelves Provided. In Case Of High Location To Avoid The Inconvenience To The Medical Staff The Pendant Can Lowered With Help Of Switch Provided On The Shelf And Operated To The Best

Some Necessary Medical Equip. Like Medical Gas Outlets And High/low Voltage Sockets. Internet Ports Are Integral Part In Medical Pendants. Also Provide IV Pole Attachments And Baskets In SS



WSP-109-B
Surgical / Anaesthesia
Pendant



Surgical Pendant Is New Addition In Field Medical Equipment In Modern Operating Theater : This Medical Pendant Has A Compact Structure By Logical Design. The Arm Console Head Can Be Rotated 340° In Low Location The Medical Staff Can Operate All The Medical Equipment In The Console Head And On The Medical Equipment Shelves Provided. In Case Of High Location To Avoid The Inconvenience To The Medical Staff The Pendant Can Lowered With Help Of Switch Provided On The Shelf And Operated To The Best

Some Necessary Medical Equip. Like Medical Gas Outlets And High/low Voltage Sockets. Internet Ports Are Integral Part In Medical Pendants. Also Provide IV Pole Attachments And Baskets In SS



WSP-110-C
Surgical / Anaesthesia
Pendant



Surgical Pendant Is New Addition In Field Medical Equipment In Modern Operating Theater : This Medical Pendant Has A Compact Structure By Logical Design. The Arm Console Head Can Be Rotated 340° The Pendant Is Manufactured In Accordance To The Ceiling Height Of The Operation Room And In Consideration With The False Ceiling To Facilitate Easy Movement Of The Pendant For The Medical Staff.

Some Necessary Medical Equip. Like Medical Gas Outlets And High/low Voltage Sockets. Internet Ports Are Integral Part In Medical Pendants. Also Provide IV Pole Attachments And Baskets In SS



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54



© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical

WSP-111-D
Single Arm
Surgical / Anaesthesia
Pendant



Single Arm Surgical Pendant Is New Addition In Field Medical Equipment In Modern Operating Theater : This Medical Pendant Has A Compact Structure By Logical Design. The Arm Console Head Can Be Rotated 340° The Pendant Is Manufactured In Accordance To The Ceiling Height Of The Operation Room And In Consideration With The False Ceiling To Facilitate Easy Movement Of The Pendant For The Medical Staff. Provided With 2 Shelves and Drawer System.

Some Necessary Medical Equip. Like Medical Gas Outlets And High/low Voltage Sockets. Internet Ports Are Integral Part In Medical Pendants. Also Provide IV Pole Attachments And Baskets In SS



WSP-112-E
Wall / Ceiling Mounted
I.C.U. Pendant



- 1) Provision For Gas Outlets
 2 Oxygen 1 Vacuum, 1 Air.
- 2) 6 Electrical Connections.
- 3) 1 RJ 45
- 4) 4 Pole IV Attachment.
- 5) Monitor Arm With Monitor Holder.
- 6) 4 Shelves For Syringe/ Infusion Pump.
- 7) Wall / Ceiling Mounted

Extra Accessories Like Gas Outlets S.S. Basket, Drawer System, Oxygen Flow Meter. Flexible Gas Tubing, Nozzles Etc. Can Also Be Provided.



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
 Email : westernsurgical@gmail.com
 Ph. +91-281-2232871 / 2243159
 Mo. +91 94280 03750/52/54

UNDER BED CRIB



WS - N13A

CRIB ON BED

WS - N13B

CRIB ON STAND

Overall Size : 35" x 27" x 40"
Tubular Structure with crib either
stainless steel or powder coated.
Supplied in KDC.



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED































Western Surgical

75

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

Accessories Hospital Furniture

			
			
			
			
	 <p>Central Locking Castors 125mm</p>	 <p>Twin Castors 125mm</p>	 <p>Heavy Castors 125mm</p>
 <p>Economy Castors 125mm</p>			
<p>WS F1 : DOCUMENT CARRIER Size : 510 x 280 x 150mm 20" x 11" x 6" S.S. 304 Grade.</p> 	 <p>ABS Expandable Food Tray</p>	 <p>Bed Side Table Top & Gas Spring</p>	 <p>Bed Crank - 48"</p>
	 <p>I.V. Attachment For Bed</p>		

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

(76)

© Copyright Reserved by : Western Surgical - Rajkot.

OUR OTHER PRODUCTS RANGE



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

OUR OTHER PRODUCTS RANGE



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.









WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

OUR OTHER PRODUCTS RANGE



BPMPR 121 Aqua (Metallic)
BPMPR 142 Wall Regular

ADULT WIGHT SCALE KRUPS EMPERIL

Bed pan Plastic M/F

Dolphine massager

Diaper,s Pediatric/Adult

Bike Air

Water Bed Jyoti / Duck Back Healing touch

Hilling pad Flamingo all range

Walking Stick & Crutches

Hearing Behind

Nebulizer - Hamant

Hot Water bottle all range

Hilling pad Flamingo all range

Walking Stick & Crutches

Hearing Behind

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.

Dt. 27-04-2024

Printing Dt. 27-04-2024



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.

Western Surgical



Western[®] Surgical

WESTERN SURGICAL EQUIPMENTS PRIVATE LIMITED

We are Available on :



GeM
Government
e Marketplace

AN ISO 9001 : 2008 & ISO 13485 : 2015 AC 2012, CE, FDA & GMP CERTIFIED COMPANY

"Happy Home", 1st Floor, 23,
Karanpara, B/h. S.T. Bus Stand,
Kanak Road, Rajkot-360001.

Ph. : +91-281-2232871/2243159,

Mo. : +91 98250 69069 - +91 98250 45997

Email : westernsurgical@gmail.com

Email : info@westernsurgical.in

* www.westernsurgical.in *

© Copyright Reserved by : Western Surgical - Rajkot.



OUR MANUFACTURING & SHOWROOM UNIT



7-Dhwani Industrial Estate,
Nr. Kothariya Railway Crossing,
B/h. Jay Ganesh Ford Service Center,
NH-8B, Gondal Road, Nr. Meera Papad,
Rajkot-360 007. Gujarat (India).